

Certificate of Conformity

VWR Catalogue Number	732-3649
Description	VWR® PCR Tubes, Classic, PCR tubes with flat caps, 0.2 mL Capacity, Natural Colour, Sterile
Lot Number	20706/2526_01
Date of Manufacture (yyyy-mm-dd)	2025-01
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2026-06-25

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 & ISO 13485:2016**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Material	Polypropylene (PP). Resins are USP Class-VI certified, FDA regulation CFR 21, and are FDA approved for food contact.
Sterilization	Products are irradiated, and dose released upon ISO 11137 recommended practices in effect at the time of validation meet a minimum requirement of 10^{-6}.
ATP Assay	Not applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination. Test limit is 1.0×10^{-6} Kunitz Units for DNase & 1.0×10^{-9} Kunitz Units for RNase.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).

Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Endotoxin levels are tested by Limulus Amoebocyte Lysate (LAL) Assay. Product is free from Endotoxin.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	This product is free of Bisphenol A.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analysis is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Room Temperature.
Shelf Life	Not applicable

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.