


## A15738 L-Arginine, 98+%

**Product Number:** A15738  
**CAS number:** 74-79-3  
**MDL number:** MFCD00002635  
**Molecular formula:** C<sub>6</sub>H<sub>14</sub>N<sub>4</sub>O<sub>2</sub>  
**Molecular weight:** 174.20



## Product Specification

**Appearance (Color):** White  
**Form:** Powder or crystalline powder  
**Assay (Non-aqueous acid-base Titration):** ≥98.0 to ≤102.0%  
**Identification (FTIR):** Conforms  
**Optical Rotation:** 26.5 ± 1.0° (C=8 in 5N HCl)

**Date of Print:** June 30, 2022  
**Version:** 9

Product Specifications are subject to amendment and may change over time.