

PRODUCT CERTIFICATE

Product: Contec Sterile NeutraKlean
Product Code: SBC103NK
Product Description: Sterile Neutral Detergent in purified water (EP) 1L Capped
Batch Number: 240402737
Manufacture Date: APR / 2024
Expiry Date: APR / 2026

ANALYSIS

Test	Specification	Results
Colour:	Colourless	Complies
Clarity:	Clear	Complies
Filtration:	Filtered to 0.2 microns	Complies
SG at 20°C:	0.995 - 1.010	1.000
pH:	6.5 to 8.0	7.58

Manufactured product via a quality system certified to ISO 9001:2015, tested in accordance with documented quality procedures and approved when required specifications are met.

IRRADIATION

Irradiation certificate number: 2173-46603A
Irradiation Dose (kGy): > 25 kGy > 28.8

We certify that the notified goods have undergone irradiation by exposure to γ (Gamma) irradiation with an exposure of not less than 25 kGy to conform to the European Pharmacopoeia.

Irradiation treatment applied was accordance with:

ISO 13485:2016 Quality Management System – Medical Devices
ISO 11137:2015 Sterilisation of Healthcare Products – Requirements for Validation & Routine Control
– Radio-sterilisation

STERILITY

Sterility test number: 0000839130
Sterility test result: No evidence of microbial growth
Test method as described in the current edition of the European Pharmacopoeia.

Name:	1: Kobey Robson	2: Nicola Hunter
Position:	1: Quality Assurance Officer	2: Quality Assurance Officer
Date:	1: 25 JUN 24	2: 27 JUN 24
Authorised Signature:	1: 	2: 

For and on behalf of Contec Inc

STERIS: Gamma Certificate Of Processing

Prepared for: CONTEC CLEANROOM UK LTD (8250)

Gamma Process Run ID 2173-46603A

<u>Product Code</u>	<u>Lot Number</u>	<u>Quantity</u>	<u>UOM</u>
1L Neutraklean DV5638	240402737	155	Case

Validation Reference Number: 5638

Processing Run Start Date: 08-May-2024 6:02 PM

Processing Run End Date: 09-May-2024 2:05 AM

Specified Dose Range (kGy):	25.0 - 45.0	Calculated Min Dose (kGy):	28.8
-----------------------------	-------------	----------------------------	------

Reference Dose Range (kGy):	30.8 - 39.3	Calculated Max Dose (kGy):	42.9
-----------------------------	-------------	----------------------------	------

PO Number: 140844

Product meets Customer specifications; zero nonconformities occurred during this irradiation run.

Gamma Process Run Approval authorized by STERIS

Date/Time: 09-May-2024 6:49 AM

Operating facilities are in compliance with applicable regulations providing services under a certified ISO 13485 quality system and processing in accordance with STERIS practices and procedures.

Processing Location

Synergy Health Sterilisation UK
Limited, a STERIS Company
Brunel Close
Drayton Fields Industrial Estate
Daventry
Northants
NN11 8RB
Phone: + 44(0) 1327 706 111


10 May 24