



Certificate of Conformity

VWR LABORATORY JACKS

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Catalog No.	Plate size (mm)	Max. static load (kg)	Max. dynamic load (kg)	Working height (mm)	Pk
245-1220	75×80	15	2	49 - 147	1
245-1221	120×140	30	5	60 - 243	1
245-1222	130×160	30	7	65 - 255	1
245-1223	200×200	30	7	70 - 260	1
245-1226	240×240	30	7	60 - 265	1
245-1227	300×300	60	12	130 - 470	1
245-1228	400×400	60	15	165 - 615	1

Material : Aluminium with stainless steel plates

Other Characteristics:

- High quality, stable and inexpensive range of lifting platforms for a variety of uses
- Plates are made from matt, brushed, non rusting stainless steel, and smooth running aluminium frame with scissor action.
- Rubber feet ensure stable support
- Ergonomically designed height adjustment knob for easy operation.
- Available with several different platform sizes.
- With M10 hole for retort rod (except 245-1220).

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **April 2017**