

Certificate of Conformity

VWR Catalogue Number	See List
Description	Translucent Euro containers, VWR®
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2022-04-18

Containers

VWR Catalogue Number	Description	Capacity [L]	W x D x H ext [mm]	W x D x H Int [mm]
216-2447	Translucent Euro container, 15 kg load capacity, closed handles	2.9	300 × 200 × 70	268 × 168 × 68
216-2448	Translucent Euro container, 15 kg load capacity, closed handles	5.1	300 × 200 × 120	268 × 168 × 118
216-2449	Translucent Euro container, 15 kg load capacity, closed handles	7.2	300 × 200 × 170	268 × 168 × 168
216-2450	Translucent Euro container, 15 kg load capacity, closed handles	9.4	300 × 200 × 220	268 × 168 × 218
216-2451	Translucent Euro container, 25 kg load capacity, open handles	20.8	400 × 300 × 220	369 × 270 × 218
216-2452	Translucent Euro container, 25 kg load capacity, closed handles	23.9	600 × 400 × 120	568 × 368 × 117
216-2453	Translucent Euro container, 25 kg load capacity, open handles	44.2	600 × 400 × 220	568 × 368 × 217
216-2454	Translucent Euro container, 25 kg load capacity, open handles	64.5	600 × 400 × 320	568 × 368 × 317

Lids

VWR Catalogue Number	Description	Dimensions [mm]
216-2455	Translucent lid 300×200, 80 kg load capacity	300 × 200 × 10
216-2456	Translucent lid 400×300, 80 kg load capacity	400 × 300 × 10
216-2457	Translucent lid 600×400, 80 kg load capacity	600 × 400 × 10

Quality System Compliance

Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Virgin polypropylene.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Store at room temperature in normal warehouse conditions.
Shelf Life	No predefined shelf life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.