

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Cleanroom Notebooks
Country of Origin	Manufactured in South Korea
Date of Issue (yyyy-mm-dd)	2025-04-25

VWR Catalogue Number	Paper Color	Sheet Size	Unit
10159-918	White	3 × 4"	Case of 40
10159-920	White	5.5 × 8.5"	Case of 20
10159-922	White	8.5 × 11"	Case of 10

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with IEST-RP-CC003.4, Appendix B, Section B.2.5 (Releasable Large Particle Test).

Case	50 sheets, 100 pages each, 20 per case	
Specifications	22 lb / 80 gsm	28 lb / 105 gsm
Basis Weight (g/m ²)	78 – 83	102 – 107
Caliper (micron)	94 – 104	110 – 120
Tear Strength MD/CD	74/75	80/81
Brightness (%)	88 – 92	88 - 92

Product Specifications

Material

Paper: Eucalyptus Pulp.

Covers: Semi-transparent, High-Density Polyethylene

Sterilization

The products are non-sterile.

ATP Assay

Not applicable.

DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Store at 48 °F to 99 °F.
Shelf Life	5 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	January 2022
2	2.0	Removed incorrect information and updated correct information	April 2025