



Certificate of Conformity

COVER CRITICAL VWR BLUE, STERILE 129-1220

CHARACTERISTICS:

- Used for equipment, instrument cover protection
- 53grams/m2 Laminated fabric
- Serged seams
- Elastic latex free
- Blue color
- Manufactured in a stringent clean environment
- Sterile (SAL-6)

PHYSICAL PROPERTIES

Material	53grams/m2 Laminated Non-Woven	
Size	Height : 137cm	
	Width: 147cm	
	Length: 66cm	
Colors	Blue	
Tensile Strength(ASTM D5034)	MD > 70N/5cm;	CD > 80N/5cm
Elongation Rate (ASTM D5034)	MD > 40%	CD > 40%
Ball Burst Strength(ASTM D3787)	>90N	
Helmke Drum Test	Level I (101 counts/ft ³ @particles ≥ 0.5um)	

PACKING:

10PCS/case

Certificate of Conformity

QUALITY ASSURANCE CRITERIA:

Protection covers are guaranteed to meet 4.0 AQL level II sampling plan per MIL-std-150E for the following defect criteria:

- | | |
|--------------------------|---|
| 1. Broken fabric | 4. Holes/tears in fabric caused by sewing machine |
| 2. Irregular PP material | 5. Physical dimensions non-com |
| 3. Dirty PP material | |

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **November 2015**