

Certificate of Conformity

Catalogue Number	See list
Description	Dispenser tips, classic, non-sterile
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2022-02-16

Catalogue Number	Item Description	Volume [mL]	UOM [pcs]
613-1000	Dispenser tips classic	0.05	100
613-1001	Dispenser tips classic	0.5	100
613-1002	Dispenser tips classic	1.25	100
613-1003	Dispenser tips classic	2.5	100
613-1004	Dispenser tips classic	5.0	100
613-1005	Dispenser tips classic	12.5	100
613-1006	Dispenser tips classic	25.0	25
613-1007	Dispenser tips classic	50.0	25
613-5886	Adapter for 25 mL tips	25	1
613-5887	Adapter for 50 mL tips	50	1
613-3811	Adapter for 50 mL tips	50	1

Quality System Compliance

Products are manufactured in **ISO 13485:2016** and **ISO 9001:2015** certified plant. Products are inspected and controlled in accordance with the current applicable product specification.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Tube: Polypropylene Adapter: Polypropylene Plunger 0.05 ml tip: high density polyethylene + pin (= stainless steel) Plunger 0.5 ml – 50.0 ml: high density polyethylene (HDPE).
Sterilization	Not applicable.
ATP Assay	Not applicable.
PCR Inhibitor	Not applicable.
DNase & RNase free	Not applicable.
Nucleic Acid Free	Not applicable.
Manufacturing Environment	Dispenser tips are manufactured under controlled conditions.
BSE/TSE	Tallow derived additives may be used in the manufacture of the raw materials. These tallow derived additives are in compliance with the guidance EMA/410/01 Rev. 3 – July 1, 2011 “Minimizing the risk of transmitting animal spongiform encephalopathy agents via human and veterinary medicinal products”, adopted by the European Commission and published in the Official Journal of the European Union March 5, 2011 (2011/C 73/01). To the best of our knowledge, these treatments ensure inactivation of any agents present in the tallow derivatives to be used in plastics materials. We therefore can state that this grade is considered safe to use with respect to BSE and TSE transmissions.
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	The raw materials used are compliant with the currently valid directive 2011/65/EU and 2002/95/EG (RoHS).
REACH Statement	The product does not contain any “SVHC” ingredients, in accordance with the currently valid regulation (EC) No 1907/2006 on REACH.

Storage Conditions Keep dry and away from sunlight.

Shelf Life 5 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.