

Certificate of Conformity

VWR Catalogue Number	See List
Description	Bulb pipettes, two mark, Class AS
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2021-12-06

Catalogue Number	Capacity [mL]	Tolerance [±mL]	Length [mm]
612-4775	1	0.007	300
612-4776	2	0.010	330
612-4777	3	0.010	350
612-4778	4	0.015	410
612-4779	5	0.015	410
612-4780	10	0.020	450
612-4781	20	0.030	520
612-4782	25	0.030	530
612-4783	50	0.050	550
612-4784	100	0.080	600

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 18001:2007** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

The product is in conformance with the requirements of the following standards and regulations

Technical standards for products: EN ISO 648

Calibration procedure: ISO 4787

Material	Borosilicate 9.1 Glass.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Room temperature.
Shelf Life	No predefined shelf life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.