




Capacity ml	Catalog number	Material	Color	Screw Closure Size, mm	Height without closure	Body Diameter	Neck ID	Autoclave able
					MM	MM	MM	
125	414004-225	LDPE	TRANSLUCENT	24-410	98	51	17.5	NA
250	414004-226	LDPE	TRANSLUCENT	24-410	131	61	17.5	NA
500	414004-227	LDPE	TRANSLUCENT	38-430	167	73	27.7	NA
1000	414004-228	LDPE	TRANSLUCENT	38-430	213	92	27.7	NA

Material Recycle Code



Description: Wash Bottle

Drawing:
Not to Scale



Release: 4/19/2022

Revision: 1