

	Document Number	Revision	Page	Date
	CSS001	A2	1 of 1	12-01-2014
PRODUCT DATA SHEET		Cell Strainers		


1: Material: Strainer: PP (Polypropylene)

Mesh: Nylon

2: Related Catalog Numbers:

VWR EU EAN. No.	VWR NA Cat. No.	Description	Color	Sterile	Packout Description
732-2757	10199-654	Cell strainers ,Fits all standard 50 ml tubes, 40µm	Blue	Yes	1/ PK, 50/Case
732-2758	10199-656	Cell strainers ,Fits all standard 50 ml tubes, 70µm	White	Yes	1/ PK, 50/Case
732-2759	10199-658	Cell strainers ,Fits all standard 50 ml tubes, 100µm	Yellow	Yes	1/ PK, 50/Case

3: Dimensions (Unit:mm)



4: Features:

- Available in 3 mesh sizes, 40 µm, 70 µm, and 100 µm
- 3 different colors: blue, white, and yellow, for easy identification
- Improved uniformity of single cell suspensions
- Made of a strong nylon mesh with evenly spaced mesh pores
- The extended lip on the strainer enables aseptic handling with forceps
- Design to fit perfectly into a 50ml centrifuge tube
- Ready-to-use, sterilized by gamma irradiation
- Individually packed in easy to open bags
- DNase- & RNase-free
- Noncytotoxic
- Non-pyrogenic

5: STERILITY:

- Sterilized by gamma irradiation SAL10⁻⁶ (ISO11137)

6: PRODUCT ISO 9000 CERTIFICATION:

- Manufactured in a class 100,000 room environment.
- Manufactured under EN ISO13485 :2003 and ISO9001:2008 quality management system.

Drawn	Approved
Date 1 Dec 2014	Date 1 Dec 2014
Name Alex Fang	Name James