

avantor[™]
delivered by **VWR**[™]


VWR.COM

AXYGEN[®]
A Corning Brand

Axygen[®] Pipet Tips

High quality pipet tips for all
your liquid handling needs





Axygen® Pipet Tips

Overview	2
Axygen MultiRack Tip System	3
Universal Filter Tips	3
Universal Tips	3
Axygen LTS Compatable Tips	5
Axygen HybridRack™ Pipet Tips	6
Axygen Traditional Pipet Tips	7
Universal Filter Tips	7
Universal Tips	8
Related Products	8


Products may not be available in all markets.

Overview

Axygen® universal pipet tips are designed to fit a wide variety of single- and multi-channel pipettors. Flexible walls and a series of internal sealing rings ensure a secure fit with less force required to load and eject.

The Axygen product offering now also includes the Axygen MultiRack System, which is our updated reload and racking system. In addition to the universal fit tips packaged in the MultiRack System, we now also offer LTS Compatible tips in the MultiRack System, and if you are worried about your carbon footprint, the Axygen HybridRack™ is a perfect option for you.

- ▶ **Automated assembly** – Ensures product is free from contamination.
- ▶ **ISO 13485 and ISO 14001** – Certified quality management systems supply stringent quality processes to ensure the highest quality finished good.
- ▶ **In-house testing** – In-house QA lab allows product to be validated for fit, form, and function prior to release.
- ▶ **Brand you can trust** – A long track record of innovation and value as one of the Corning trusted brands.
- ▶ **Broad product offering** – Corning offers many different styles of pipet tips, including racked, filtered, reloads, stack racks, low retention, wide bore, and more. With so many options, you are sure to find a tip that works for you.



Packaging Solutions




Bulk



Reload



Racked



Tip Station



HybridRack

Corning products are included in Avantor's Environmentally Preferable Products (EPP) program, designed to provide greater transparency into product sustainability attributes. This allows you to make more informed and sustainable purchasing decisions.

Axygen® MultiRack Tip System



Engineered to fit most pipettors on the market, the Axxygen MultiRack Tip System is a complete pipet tip system offering in various packaging styles to meet every lab's needs. The MultiRack Tip System offers standard and extended-length sizes in volumes from 0.1 µL to 1250 µL. Racked and reload packaging styles feature preloaded tips in easy-load, color-coded deck inserts. With a small benchtop footprint and reusable rack, this system minimizes environmental impact and maximizes lab space. Axxygen MultiRack Tip System products are made in the USA and meet ISO standards.

Axygen MultiRack Universal Filter Pipet Tips

- ▶ Racked
- ▶ Sterile
- ▶ RNase-/DNase-free
- ▶ Human DNA-free
- ▶ PCR inhibitor-free
- ▶ Non-pyrogenic
- ▶ USP Class VI polypropylene



Standard Surface VWR Cat. No. Corning Cat. No.	Maxymum Recovery Surface VWR Cat. No. Corning Cat. No.	Volume	Description	Quantity
76320-724 MRF-10XT-R-S	76320-750 MRF-10XT-L-R-S	10 µL	Extended length, graduated, fine point	96/rack, 10 racks/pack, 5 packs/case
76320-726 MRF-20-R-S	76320-752 MRF-20-L-R-S	20 µL	Graduated, fine point	96/rack, 10 racks/pack, 5 packs/case
76320-728 MRF-100-R-S	76320-754 MRF-100-L-R-S	100 µL	Graduated, fine point	96/rack, 10 racks/pack, 5 packs/case
76320-730 MRF-200NX-R-S	76320-736 MRF-200NX-L-R-S	200 µL	Graduated, fine point	96/rack, 10 racks/pack, 5 packs/case
76320-732 MRF-200XT-R-S	76320-738 MRF-200XT-L-R-S	200 µL	Extended length, graduated, fine point	96/rack, 8 racks/pack, 5 packs/case
76320-734 MRF-300XT-R-S	76320-740 MRF-300XT-L-R-S	300 µL	Extended length, graduated, fine point	96/rack, 8 racks/pack, 5 packs/case
76320-748 MRF-1000XT-R-S	76320-742 MRF-1000XT-L-R-S	1250 µL	Extended length, graduated, fine point	96/rack, 8 racks/pack, 5 packs/case

Axygen MultiRack Universal Pipet Tips

10 µL Extended Length, Graduated, Fine Point



Standard Surface VWR Cat. No. Corning Cat. No.	Maxymum Recovery Surface VWR Cat. No. Corning Cat. No.	Packaging	Sterile	Quantity
76320-720 MR-10XT	-	Bulk	No	1,000/pack, 10 packs/case
76320-714 MR-10XT-RL	-	Reload	No	96/deck, 10 decks/pack, 4 packs/case
76320-676 MR-10XT-R	-	Racked	No	96/rack, 10 racks/pack, 5 packs/case
76320-686 MR-10XT-R-S	76320-696 MR-10XT-L-R-S	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
76320-706 MR-10XT-TS	-	Tip station	No	192/layer, 5 layers/pack, 5 packs/case
76320-708 MR-10XT-TS-S	-	Tip station	Yes	192/layer, 5 layers/pack, 5 packs/case

Axygen® MultiRack Universal Pipet Tips (continued)

200 µL, Graduated, Fine Point

Standard Surface VWR Cat. No. Corning Cat. No.	Maxymum Recovery® Surface VWR Cat. No. Corning Cat. No.	Packaging	Sterile	Quantity
76320-672 MR-200	-	Bulk	No	1,000/pack, 10 packs/case
76320-716 MR-200-RL	-	Reload	No	96/deck, 10 decks/pack, 4 packs/case
76320-678 MR-200-R	-	Racked	No	96/rack, 10 racks/pack, 5 packs/case
76320-688 MR-200-R-S	76320-698 MR-200-L-R-S	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
76320-710 MR-200-TS	-	Tip station	No	192/layer, 5 layers/pack, 5 packs/case
76320-712 MR-200-TS-S	-	Tip station	Yes	192/layer, 5 layers/pack, 5 packs/case

300 µL, Graduated, Fine Point

76320-718 MR-300-RL	-	Reload	No	96/deck, 10 decks/pack, 4 packs/case
76320-680 MR-300-R	-	Racked	No	96/rack, 10 racks/pack, 5 packs/case
76320-690 MR-300-R-S	76320-700 MR-300-L-R-S	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case

300 µL Graduated, Extended Length, Fine Point

76320-682 MR-300XT-R	-	Racked	No	96/rack, 8 racks/pack, 5 packs/case
76320-692 MR-300XT-R-S	76320-702 MR-300XT-L-R-S	Racked	Yes	96/rack, 8 racks/pack, 5 packs/case

1250 µL Graduated, Extended Length, Fine Point

76320-674 MR-1000XT	-	Bulk	No	768/pack, 10 packs/case
76320-722 MR-1000XT-RL	-	Reload	No	96/deck, 8 decks/pack, 4 packs/case
76320-684 MR-1000XT-R	-	Racked	No	96/rack, 8 racks/pack, 5 packs/case
76320-694 MR-1000XT-R-S	76320-704 MR-1000XT-L-R-S	Racked	Yes	96/rack, 8 racks/pack, 5 packs/case

Axygen® LTS Compatible Tips



Axygen LTS Compatible Tips are ergonomically designed to reduce discomfort and optimized to work with Rainin LiteTouch pipettors, all with the quality you have come to expect from Corning. By reducing the force required to seal and eject, Axygen LTS Compatible tips can help you avoid using excess force to obtain an effective seal between the pipettor and pipet tip. Axygen LTS Compatible tips are available in your favorite packaging formats including Axygen HybridRack™, a sustainable packaging choice. All LTS tips are composed of USP Class VI polypropylene and manufactured in an ISO 13485 and ISO 14001 United States based facility.

- ▶ RNase-/DNase-free
- ▶ Human gDNA-free
- ▶ PCR inhibition-free
- ▶ Nonpyrogenic
- ▶ USP Class VI polypropylene

20 µL LTS Compatible Extended Length, Graduated, Fine Point



VWR Cat. No.	Corning Cat. No.	Packaging	Sterile	Quantity
77395-566	LTF-20-R-S	Racked, filtered	Yes	96/rack, 10 racks/pack, 5 packs/case
77395-598	LTFER-20-R-S	HybridRack, filtered	Yes	96/rack, 10 racks/pack, 5 packs/case
77395-572	LT-20-R-S	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
77395-578	LT-20-R	Racked	No	96/rack, 10 racks/pack, 5 packs/case
77395-584	LT-20-RL	Reload stack	No	96/deck, 10 decks/pack, 4 packs/case
77395-590	LT-20-TS	Tip station	No	192/layer, 5 layers/pack, 5 packs/case
77395-592	LT-20-TS-S	Tip station	Yes	192/layer, 5 layers/pack, 5 packs/case

200 µL LTS Compatible, Graduated, Fine Point



VWR Cat. No.	Corning Cat. No.	Packaging	Sterile	Quantity
77395-568	LTF-200-R-S	Racked, filtered	Yes	96/rack, 10 racks/pack, 5 packs/case
77395-600	LTFER-200-R-S	HybridRack, filtered	Yes	96/rack, 10 racks/pack, 5 packs/case
77395-574	LT-200-R-S	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
77395-580	LT-200-R	Racked	No	96/rack, 10 racks/pack, 5 packs/case
77395-586	LT-200-RL	Reload stack	No	96/deck, 10 decks/pack, 4 packs/case
77395-594	LT-200-TS	Tip station	No	192/layer, 5 layers/pack, 5 packs/case
77395-596	LT-200-TS-S	Tip station	Yes	192/layer, 5 layers/pack, 5 packs/case

1000 µL LTS Compatible, Graduated, Fine Point



VWR Cat. No.	Corning Cat. No.	Packaging	Sterile	Quantity
77395-570	LTF-1000-R-S	Racked, filtered	Yes	96/rack, 8 racks/pack, 5 packs/case
77395-602	LTFER-1000-R-S	HybridRack, filtered	Yes	96/rack, 10 racks/pack, 5 packs/case
77395-576	LT-1000-R-S	Racked	Yes	96/rack, 8 racks/pack, 5 packs/case
77395-582	LT-1000-R	Racked	No	96/rack, 8 racks/pack, 5 packs/case
77395-588	LT-1000-RL	Reload stack	No	96/deck, 8 decks/pack, 4 packs/case

Axygen® HybridRack™ Pipet Tips



Born from a desire to reduce the amount of waste happening in laboratories, Axxygen HybridRack pipet tips are an example of our innovative spirit leading to a more sustainable future. By replacing the traditional plastic rack and lid normally found in racked pipet tip products with a paperboard rack and lid, we can reduce the amount of plastic per case by 70%, or approximately two pounds per case of 200 μ L tips. The HybridRack is strong enough to support multi-channel pipetting, and as soon as the rack is empty, it can be recycled using any curbside recycling program.

- ▶ Sterile
- ▶ RNase-/DNase-free
- ▶ Human DNA-free
- ▶ PCR inhibitor-free
- ▶ Non-pyrogenic
- ▶ USP Class VI polypropylene

10 μ L Extended Length, Graduated, Fine Point

Standard Surface VWR Cat. No. Corning Cat. No.	Maxymum Recovery® Surface VWR Cat. No. Corning Cat. No.	Filtered	Sterile	Quantity
77269-858 ERF-10XT-R-S	77269-852 ERF-10XT-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77269-846 ER-10XT-R-S	-	No	Yes	96/rack, 10 racks/pack, 5 packs/case

20 μ L, 100 μ L, and 200 μ L Graduated, Fine Point

77269-866 ERF-20-R-S	77269-864 ERF-20-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77269-868 ERF-100-R-S	77269-860 ERF-100-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77269-876 ER-200-R-S	-	No	Yes	96/rack, 10 racks/pack, 5 packs/case

200 μ L and 300 μ L Graduated, Fine Point

77269-870 ERF-200NX-R-S	77269-862 ERF-200NX-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77269-878 ER-300NX-R-S	-	No	Yes	96/rack, 10 racks/pack, 5 packs/case

300 μ L Graduated, Extended Length, Fine Point

77269-872 ERF-300XT-R-S	77269-854 ERF-300XT-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77269-848 ER-300XT-R-S	-	No	Yes	96/rack, 10 racks/pack, 5 packs/case


1250 μ L Graduated, Extended Length, Fine Point

77269-874 ERF-1000XT-R-S	77269-856 ERF-1000XT-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77269-850 ER-1000XT-R-S	-	No	Yes	96/rack, 10 racks/pack, 5 packs/case

To learn more about the Corning HybridRack tips, visit vwr.com/corning.



Axygen® Traditional Pipet Tips



Known around the world for industry leading quality, accuracy, and precision, Axxygen Traditional pipet tips have been the preferred liquid handling solution for laboratories around the world for the last 20+ years. Axxygen Traditional pipet tips are designed to fit a wide variety of single- and multi-channel pipettors. Flexible walls and a series of internal sealing rings ensure a secure fit with less force required to load and eject. Featuring tips with extended length, wide bore designs for viscous solutions, and our trademarked Maxymum Recovery surface, Axxygen Traditional pipet tips are sure to meet and exceed your liquid handling expectations.

- ▶ RNase-/DNase-free
- ▶ Human gDNA-free
- ▶ Nonpyrogenic
- ▶ USP Class VI polypropylene

Axygen Traditional Universal Filter Pipet Tips

- ▶ Racked
- ▶ Filtered
- ▶ Sterile



Standard Surface VWR Cat. No. Corning Cat. No.	Maxymum Recovery® Surface VWR Cat. No. Corning Cat. No.	Volume (μL)	Description	Quantity
10011-360 TF-300-R-S	22234-018 TF-300-L-R-S	0.5 - 10	Graduated, fine point	96/rack, 10 racks/pack, 5 packs/case
10011-364 TF-400-R-S	22234-020 TF-400-L-R-S	0.5 - 10	Extended length, graduated, beveled	96/rack, 10 racks/pack, 5 packs/case
47749-872 TF-420-R-S	47749-870 TF-420-L-R-S	0.5 - 20	Extended length, graduated, beveled	96/rack, 10 racks/pack, 5 packs/case
47745-492 TXLF-10-R-S	47749-868 TXLF-10-L-R-S	0.5 - 10	Extended length, fine point	96/rack, 10 racks/pack, 5 packs/case
10011-358 TF-20-R-S	22234-008 TF-20-L-R-S	0.5 - 20	Beveled	96/rack, 10 racks/pack, 5 packs/case
10011-370 TF-50-R-S	22234-010 TF-50-L-R-S	0.5 - 50	Graduated, beveled	96/rack, 10 racks/pack, 5 packs/case
10011-352 TF-100-R-S	22234-012 TF-100-L-R-S	1 - 100	Graduated, beveled	96/rack, 10 racks/pack, 5 packs/case
10011-356 TF-200-R-S	22234-016 TF-200-L-R-S	1 - 200	Graduated, beveled	96/rack, 10 racks/pack, 5 packs/case
80082-898 TF-205-WB-R-S	80082-902 TF-205-WB-L-R-S	1 - 200	Wide bore	96/rack, 10 racks/pack, 5 packs/case
10011-362 TF-350-R-S	22234-132 TF-350-L-R-S	1 - 300	Fine point	96/rack, 10 racks/pack, 5 packs/case
10011-350 TF-1000-R-S	22234-024 TF-1000-L-R-S	100 - 1000	Beveled	100/rack, 10 racks/pack, 5 packs/case
80082-900 TF-1005-WB-R-S	80082-904 TF-1005-WB-L-R-S	1 - 1000	Wide bore	100/rack, 10 racks/pack, 5 packs/case

Axygen® Traditional Universal Pipet Tips

10 µL Graduated, Fine Point



Standard Surface VWR Cat. No. Corning Cat. No.	Maximum Recovery® Surface VWR Cat. No. Corning Cat. No.	Volume (µL)	Packaging	Sterile	Quantity
10011-928 T-300	22234-092 T-300-L	0.5 - 10	Bulk	No	1,000/pack, 20 packs/case
10011-336 T-300-R	22234-094 T-300-L-R	0.5 - 10	Racked	No	96/rack, 10 racks/pack, 5 packs/case
10011-338 T-300-R-S	22234-096 T-300-L-R-S	0.5 - 10	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
10011-874 T-300-STK	10011-872 T-300-L-STK	0.5 - 10	Stack rack	No	96/layer, 5 layers/stack, 10 stacks/case
47734-724 T-300-STK-S	47734-740 T-300-L-STK-S	0.5 - 10	Stack rack	Yes	96/layer, 5 layers/stack, 10 stacks/case
89029-921 RFL-300-C	-	0.5 - 10	Reload	No	96 tips/deck, 10 decks/pack, 5 packs/case

10 µL Extended Length, Graduated, Beveled



10011-930 T-400	22234-098 T-400-L	0.5 - 20	Bulk	No	1,000/pack, 20 packs/case
10011-344 T-400-R	22234-100 T-400-L-R	0.5 - 20	Racked	No	96/rack, 10 racks/pack, 5 packs/case
10011-346 T-400-R-S	22234-102 T-400-L-R-S	0.5 - 20	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case

10 µL Extended Length, Fine Point



47745-484 TXL-10	12000-208 TXL-10-L	0.5 - 20	Bulk	No	1,000/pack, 20 packs/case
47745-486 TXL-10-R	12000-210 TXL-10-L-R	0.5 - 20	Racked	No	96/rack, 10 racks/pack, 5 packs/case
47745-488 TXL-10-R-S	12000-212 TXL-10-L-R-S	0.5 - 20	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case

200 µL Beveled



10011-918 T-200-C	22234-062 T-200-C-L	1 - 200	Bulk	No	1,000/pack, 20 packs/case
10011-316 T-200-C-R	22234-064 T-200-C-L-R	1 - 200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
10011-318 T-200-C-R-S	22234-066 T-200-C-L-R-S	1 - 200	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
10011-868 T-200-C-STK	22234-068 T-200-C-L-STK	1 - 200	Stack rack	No	96/layer, 5 layers/stack, 10 stacks/case

200 µL Beveled, Yellow



10011-920 T-200-Y	-	1 - 200	Bulk	No	1,000/pack, 20 packs/case
10011-320 T-200-Y-R	-	1 - 200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
10011-322 T-200-Y-R-S	-	1 - 200	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
10011-870 T-200-Y-STK	-	1 - 200	Stack rack	No	96/layer, 5 layers/stack, 10 stacks/case

Axygen® Traditional Universal Pipet Tips (continued)

200 µL Reference-compatible, Beveled, Graduated



Standard Surface VWR Cat. No. Corning Cat. No.	Maxymum Recovery® Surface VWR Cat. No. Corning Cat. No.	Volume (µL)	Packaging	Sterile	Quantity
10011-932 TR-222-C	22234-070 TR-222-C-L	1 - 200	Bulk	No	1,000/pack, 20 packs/case
10011-372 TR-222-C-R	22234-072 TR-222-C-L-R	1 - 200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
10011-374 TR-222-C-R-S	22234-074 TR-222-C-L-R-S	1 - 200	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
10011-894 TR-222-C-STK	22234-076 TR-222-C-L-STK	1 - 200	Stack rack	No	96/layer, 5 layers/stack, 10 stacks/case
95042-896 TR-222-C-STK-S	47734-742 TR-222-C-L-STK-S	1 - 200	Stack rack	Yes	96/layer, 5 layers/stack, 10 stacks/case

200 µL Reference-compatible, Beveled, Graduated, Yellow



10011-934 TR-222-Y	-	1 - 200	Bulk	No	1000/pack, 20 packs/case
10011-376 TR-222-Y-R	-	1 - 200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
10011-378 TR-222-Y-R-S	-	1 - 200	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
47734-722 TR-222-Y-STK	-	1 - 200	Stack rack	No	96/layer, 5 layers/stack, 10 stacks/case
89029-919 RFL-222-C	-	1 - 200	Reload	No	96 tips/deck, 10 decks/pack, 5 packs/case
89029-922 RFL-R	-	-	Empty Rack for 89029-921 and 89029-919	No	10 racks/pack, 5 packs/case

200 µL Wide Bore



10011-922 T-205-WB-C	-	1 - 200	Bulk	No	1,000/pack, 20 packs/case
10011-324 T-205-WB-C-R	-	1 - 200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
10011-326 T-205-WB-C-R-S	-	1 - 200	Racked	Sterile	96/rack, 10 racks/pack, 5 packs/case

300 µL Fine Point



10011-876 T-350-C	10011-011 T-350-C-L	1 - 300	Bulk	No	1,000/pack, 10 packs/case
10011-340 T-350-C-R	22234-130 T-350-C-L-R	1 - 300	Racked	No	96/rack, 10 racks/pack, 5 packs/case
10011-342 T-350-C-R-S	22234-128 T-350-C-L-R-S	1 - 300	Racked	Sterile	96/rack, 10 racks/pack, 5 packs/case

1000 µL Beveled, Clear



10011-306 T-1000-C	22234-086 T-1000-C-L	100 - 1000	Bulk	No	1,000/pack, 5 packs/case
10011-308 T-1000-C-R	22234-088 T-1000-C-L-R	100 - 1000	Racked	No	100/rack, 10 racks/pack, 5 packs/case
10011-310 T-1000-C-R-S	22234-090 T-1000-C-L-R-S	100 - 1000	Racked	Yes	100/rack, 10 racks/pack, 5 packs/case

Axygen® Traditional Universal Pipet Tips (continued)

1000 µL Beveled, Blue



Standard Surface VWR Cat. No. Corning Cat. No.	Maxymum Recovery® Surface VWR Cat. No. Corning Cat. No.	Volume (µL)	Packaging	Sterile	Quantity
10011-300 T-1000-B	-	100 - 1000	Bulk	No	1,000/pack, 5 packs/case
10011-302 T-1000-B-R	-	100 - 1000	Racked	No	100/rack, 10 racks/pack, 5 packs/case
10011-304 T-1000-B-R-S	-	100 - 1000	Racked	Yes	100/rack, 10 racks/pack, 5 packs/case
89080-484 RFL-1200-C	-	1 - 1200	Reload	No	96 tips/deck, 8 decks/pack, 6 packs/case
89080-486 RFL-R-1200	-	-	Empty Rack for 89080-484	No	10 racks/pack, 3 packs/case


1000 µL Wide Bore



10011-916 T-1005-WB-C	-	100 - 1000	Bulk	No	1,000/pack, 5 packs/case
10011-312 T-1005-WB-C-R	47734-728 T-1005-WB-C-L-R	100 - 1000	Racked	No	100/rack, 10 racks/pack, 5 packs/case
10011-314 T-1005-WB-C-R-S	-	100 - 1000	Racked	Yes	100/rack, 10 racks/pack, 5 packs/case

Related Products

Axygen SnapLock Microcentrifuge Tubes



VWR Cat. No.	Corning Cat. No.	Description	Quantity
10011-680	MCT-060-C	0.6 mL Snaplock microcentrifuge tubes, MaxyClear	1,000/pack, 10 packs/case
10011-682	MCT-060-C-S	0.6 mL Snaplock microcentrifuge tubes, MaxyClear, sterile	100/bag, 5 bags/pack, 10 packs/case
10011-700	MCT-150-C	1.5 mL Snaplock microcentrifuge tubes, MaxyClear	5,000/pack, 10 packs/case
10011-702	MCT-150-C-S	1.5 mL Snaplock microcentrifuge tubes, MaxyClear, sterile	50/bag, 5 bags/pack, 10 packs/case
10011-722	MCT-175-C	1.7 mL Snaplock microcentrifuge tubes, MaxyClear	500/pack, 10 packs/case
10011-724	MCT-175-C-S	1.7 mL Snaplock microcentrifuge tubes, MaxyClear, sterile	250/pack, 10 packs/case
10011-742	MCT-200-C	2.0 mL Snaplock microcentrifuge tubes, MaxyClear	5,000/pack, 10 packs/case
10011-744	MCT-200-C-S	2.0 mL Snaplock microcentrifuge tubes, MaxyClear, sterile	50/bag, 5 bags/pack, 10 packs/case
89425-665	MCT-500-C	5.0 mL Snaplock microcentrifuge tubes, MaxyClear	250/pack, 5 packs/case
89425-667	MCT-500-C-S	5.0 mL Snaplock microcentrifuge tubes, MaxyClear, sterile	250/pack, 5 packs/case

NOTE: Other options are available, including full cases of different colors, assorted colors in the same case, capless tubes, and tubes with Axxygen Maxymum Recovery surface.

For more related products, visit vwr.com/corning.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

For additional product or technical information, visit vwr.com/corning, call 1.800.932.5000, or contact your Avantor representative.



VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2024 Avantor, Inc. All rights reserved.