



## PRODUCT DATA SHEET

# DISPOSABLE POT FOR ANATOMICAL SPECIMEN

### Description

Disposable container in PS polystyrene with sealing closure, for the possible storage of anatomical specimen.

### Intended use and mark:

Product designed for the conservation (for example in formalin) of biopsy or autopsy samples from the human body to be analyzed, in order to provide information on the physiological or pathological state of a patient, according to the European Directive 98/79 / EC on In Vitro Diagnostic medical devices (IVDD) Annex III, for which it bears the marking .

### Features

Made of shockproof PS polystyrene with high mechanical resistance, it is light, sturdy and stackable. The cap is made of PE polyethylene and its special closure guarantees, in normal use, an excellent seal against liquids and vapors. It is available in a serigraphy version with the symbols relating to biological risk and in a neutral, white version.

For capacities of 2.5 and 5 liters it is possible to mount a practical "double" handle, which facilitates the transport of the container.

### Technical and packaging data

| Code   |                                                                                                | Volume<br>(ml) | Dimension<br>Ø x H (mm) | Packaging |
|--------|------------------------------------------------------------------------------------------------|----------------|-------------------------|-----------|
| White  | Serigraphy  |                |                         |           |
| 190301 | 190300                                                                                         | 250            | 125 x 54                | 240 pz.   |
| 190501 | 190500                                                                                         | 500            | 125 x 70                | 180 pz.   |
| 191001 | 191000                                                                                         | 1000           | 125 x 120               | 180 pz.   |
| 192501 | 192500                                                                                         | 2500           | 200 x 130               | 52 pz.    |
| 195001 | 195000                                                                                         | 5000           | 200 x 215               | 40 pz.    |



| Handle | Compatibility                     |
|--------|-----------------------------------|
| 198313 | 192501 - 192500 - 195001 - 195000 |



### Environmental information:

This product is manufactured using raw materials which are environmentally friendly and fully recyclable. If necessary, dispose of this product in accordance with applicable local waste treatment Regulations

*All dimensions and volumes indicated are nominal. Technical specifications are subject to change without prior notice.*