

Tubes for single-channel pumpheads (incl. RFID transponder)

Order number	Article designation	ID mm	wall thickness mm	Lenght m	PU per piece
952 12 016*	Tygon® E-LFL	6.4	1.6	7.5	1
952 13 010	PharMed® BPT	0.8	1.6	7.5	1
952 13 011	PharMed® BPT	1.6	1.6	7.5	1
952 13 012	PharMed® BPT	2.4	1.6	7.5	1
952 13 013	PharMed® BPT	3.2	1.6	7.5	1
952 13 015	PharMed® BPT	4.8	1.6	7.5	1
952 13 016	PharMed® BPT	6.4	1.6	7.5	1
952 13 017	PharMed® BPT	8	1.6	7.5	1
952 14 011	Versilon™ 2001	1.6	1.6	15	1
952 14 013	Versilon™ 2001	3.2	1.6	15	1
952 14 015	Versilon™ 2001	4.8	1.6	15	1
952 14 016	Versilon™ 2001	6.4	1.6	15	1
952 14 017	Versilon™ 2001	8	1.6	15	1
952 14 018	Versilon™ 2001	9.5	1.6	15	1
952 15 013	Norprene® Chemical	3.2	1.6	15	1
952 15 015	Norprene® Chemical	4.8	1.6	15	1
952 15 016	Norprene® Chemical	6.4	1.6	15	1
952 15 018	Norprene® Chemical	9.5	1.6	15	1
952 17 011	Tygon® 3350	1.6	1.6	15	1
952 17 012	Tygon® 3350	2.4	1.6	15	1
952 17 013	Tygon® 3350	3.2	1.6	15	1
952 17 015	Tygon® 3350	4.8	1.6	15	1
952 17 016	Tygon® 3350	6.4	1.6	15	1
952 17 017	Tygon® 3350	8	1.6	15	1
952 17 018	Tygon® 3350	9.5	1.6	15	1
952 17 019	Tygon® 3350	11.2	1.6	15	1
952 17 035	Tygon® 3350	4.8	2.4	15	1
952 17 036	Tygon® 3350	6.4	2.4	15	1
952 17 037	Tygon® 3350	8	2.4	15	1
952 17 038	Tygon® 3350	9.5	2.4	15	1
952 18 011	Silikon Peroxide	1.6	1.6	15	1
952 18 013	Silikon Peroxide	3.2	1.6	15	1
952 18 015	Silikon Peroxide	4.8	1.6	15	1
952 18 016	Silikon Peroxide	6.4	1.6	15	1
952 18 017	Silikon Peroxide	8	1.6	15	1
952 18 035	Silikon Peroxide	4.8	2.4	15	1
952 18 036	Silikon Peroxide	6.4	2.4	15	1
952 19 013	Tygon® F-4040-A	3.2	1.6	15	1
952 19 015	Tygon® F-4040-A	4.8	1.6	15	1
952 19 016	Tygon® F-4040-A	6.4	1.6	15	1
952 19 017	Tygon® F-4040-A	8	1.6	15	1
952 20 013	Fluran® F-5500-A	3.2	1.6	15	1
952 20 015	Fluran® F-5500-A	4.8	1.6	15	1
952 20 016	Fluran® F-5500-A	6.4	1.6	15	1
952 20 017	Fluran® F-5500-A	8	1.6	15	1

Tubes for single-channel pumpheads (incl. RFID transponder)

Order number	Article designation	ID mm	wall thickness mm	Lenght m	PU per piece
952 21 011	Norprene® A-60-G	1.6	1.6	15	1
952 21 013	Norprene® A-60-G	3.2	1.6	15	1
952 21 015	Norprene® A-60-G	4.8	1.6	15	1
952 21 016	Norprene® A-60-G	6.4	1.6	15	1
952 21 017	Norprene® A-60-G	8	1.6	15	1
952 21 018	Norprene® A-60-G	9.5	1.6	15	1
952 21 019	Norprene® A-60-G	11.1	1.6	15	1
952 21 035	Norprene® A-60-G	4.8	2.4	15	1
952 21 036	Norprene® A-60-G	6.4	2.4	15	1
952 21 037	Norprene® A-60-G	8	2.4	15	1
952 21 038	Norprene® A-60-G	9.5	2.4	15	1
952 22 011	Tygon® LMT-55	1.6	1.6	15	1
952 22 012	Tygon® LMT-55	2.4	1.6	15	1
952 22 013	Tygon® LMT-55	3.2	1.6	15	1
952 22 014	Tygon® LMT-55	4	1.6	15	1
952 22 015	Tygon® LMT-55	4.8	1.6	15	1
952 22 016	Tygon® LMT-55	6.4	1.6	15	1
952 22 017	Tygon® LMT-55	8	1.6	15	1
952 22 018	Tygon® LMT-55	9.6	1.6	15	1
952 22 019	Tygon® LMT-55	11.2	1.6	15	1
952 22 035	Tygon® LMT-55	4.8	2.4	15	1
952 22 036	Tygon® LMT-55	6.4	2.4	15	1
952 22 037	Tygon® LMT-55	8	2.4	15	1
952 22 038	Tygon® LMT-55	9.6	2.4	15	1
952 23 011	Tygon® S3 E-LFL	1.6	1.6	7.5	1
952 23 013	Tygon® S3 E-LFL	3.2	1.6	7.5	1
952 23 015	Tygon® S3 E-LFL	4.8	1.6	7.5	1
952 23 016**	Tygon® S3 E-LFL	6.4	1.6	7.5	1
952 23 017	Tygon® S3 E-LFL	8	1.6	7.5	1

Stopper tubes for multi-channel pumpheads (incl. RFID transponder)

Order number	Article designation	ID mm	Color Code	wall thickness mm	Lenght m	PU per piece
952 13 703	PharMed® BPT 2 Stopper 104	0.25	orange-blue	0.91	381	6
952 13 710	PharMed® BPT 2 Stopper 104	0.89	orange-orange	0.86	381	6
952 13 717	PharMed® BPT 2 Stopper 104	1.42	yellow-yellow	0.86	381	6
952 13 722	PharMed® BPT 2 Stopper 104	2.06	violet-violet	0.86	381	6
952 13 725	PharMed® BPT 2 Stopper 104	2.79	violet-white	0.86	381	6
952 14 704	Versilon™ 2001 2 Stopper 104	0.38	orange-grün	0.91	381	6
952 14 708	Versilon™ 2001 2 Stopper 104	0.64	orange-white	0.91	381	6
952 14 712	Versilon™ 2001 2 Stopper 104	1.02	white-white	0.86	381	6
952 14 718	Versilon™ 2001 2 Stopper 104	1.52	yellow-blue	0.86	381	6
952 14 722	Versilon™ 2001 2 Stopper 104	2.06	violet-violet	0.86	381	6
952 14 725	Versilon™ 2001 2 Stopper 104	2.79	violet-white	0.86	381	6
952 17 710	Tygon® 3350 2 Stopper 104	0.89	orange-orange	0.86	381	6
952 17 717	Tygon® 3350 2 Stopper 104	1.42	yellow-yellow	0.86	381	6
952 17 722	Tygon® 3350 2 Stopper 104	2.06	violet-violet	0.86	381	6
952 17 725	Tygon® 3350 2 Stopper 104	2.79	violet-white	0.86	381	6
952 19 703	Tygon® F4040A 2 Stopper 104	0.25	orange-blue	0.91	381	12
952 19 710	Tygon® F4040A 2 Stopper 104	0.89	orange-orange	0.86	381	12
952 19 717	Tygon® F4040A 2 Stopper 104	1.42	yellow-yellow	0.86	381	12
952 19 722	Tygon® F4040A 2 Stopper 104	2.06	violet-violet	0.86	381	12
952 19 725	Tygon® F4040A 2 Stopper 104	2.79	violet-white	0.86	381	12
952 20 710	Fluran® F5500A 2 Stopper 104	0.89	orange-orange	0.86	381	12
952 20 717	Fluran® F5500A 2 Stopper 104	1.42	yellow-yellow	0.86	381	12
952 20 722	Fluran® F5500A 2 Stopper 104	2.06	violet-violet	0.86	381	12
952 20 725	Fluran® F5500A 2 Stopper 104	2.79	violet-white	0.86	381	12
952 22 702	Tygon® LMT-55 2 Stopper 104	0.19	orange-red	0.91	381	12
952 22 703	Tygon® LMT-55 2 Stopper 104	0.25	orange-blue	0.91	381	12
952 22 710	Tygon® LMT-55 2 Stopper 104	0.89	orange-orange	0.86	381	12
952 22 717	Tygon® LMT-55 2 Stopper 104	1.42	yellow-yellow	0.86	381	12
952 22 722	Tygon® LMT-55 2 Stopper 104	2.06	violet-violet	0.86	381	12
952 22 725	Tygon® LMT-55 2 Stopper 104	2.79	violet-white	0.86	381	12
952 23 703	Tygon® S3 E-LFL 2 Stopper 104	0.25	orange-blue	0.91	381	12
952 23 706	Tygon® S3 E-LFL 2 Stopper 104	0.51	orange-yellow	0.91	381	12
952 23 708	Tygon® S3 E-LFL 2 Stopper 104	0.64	orange-white	0.91	381	12
952 23 710	Tygon® S3 E-LFL 2 Stopper 104	0.89	orange-orange	0.86	381	12
952 23 712	Tygon® S3 E-LFL 2 Stopper 104	1.02	white-white	0.86	381	12
952 23 718	Tygon® S3 E-LFL 2 Stopper 104	1.52	yellow-blue	0.86	381	12
952 23 722	Tygon® S3 E-LFL 2 Stopper 104	2.06	violet-violet	0.86	381	12
952 23 725	Tygon® S3 E-LFL 2 Stopper 104	2.79	violet-white	0.86	381	12

Extension tubes for stoppertubes (incl. RFID transponder)

Order number	Article designation	ID mm	wall thickness mm	Lenght m	PU per piece
952 13 503	PharMed® BPT extension tube	0.25	0.91	3	1
952 13 506	PharMed® BPT extension tube	0.51	0.91	3	1
952 13 509	PharMed® BPT extension tube	0.76	0.91	3	1
952 13 510	PharMed® BPT extension tube	0.89	0.86	3	1
952 13 512	PharMed® BPT extension tube	1.02	0.86	3	1
952 13 516	PharMed® BPT extension tube	1.30	0.86	3	1
952 13 517	PharMed® BPT extension tube	1.42	0.86	3	1
952 13 522	PharMed® BPT extension tube	2.06	0.86	3	1
952 13 525	PharMed® BPT extension tube	2.79	0.86	3	1

952 14 504	Versilon™ 2001 extension tube	0.38	0.91	10	1
952 14 512	Versilon™ 2001 extension tube	1.02	0.86	10	1
952 14 518	Versilon™ 2001 extension tube	1.52	0.86	10	1
952 14 522	Versilon™ 2001 extension tube	2.06	0.86	10	1
952 14 525	Versilon™ 2001 extension tube	2.79	0.86	10	1

952 17 502	Tygon® 3350 extension tube	0.19	0.91	10	1
952 17 503	Tygon® 3350 extension tube	0.25	0.91	10	1
952 17 510	Tygon® 3350 extension tube	0.89	0.86	10	1
952 17 517	Tygon® 3350 extension tube	1.42	0.86	10	1
952 17 522	Tygon® 3350 extension tube	2.06	0.86	10	1
952 17 525	Tygon® 3350 extension tube	2.79	0.86	10	1

952 19 503	Tygon® F-4040-A extension tube	0.25	0.91	3	1
952 19 510	Tygon® F-4040-A extension tube	0.89	0.86	3	1
952 19 517	Tygon® F-4040-A extension tube	1.42	0.86	3	1
952 19 522	Tygon® F-4040-A extension tube	2.06	0.86	3	1
952 19 525	Tygon® F-4040-A extension tube	2.79	0.86	3	1

952 20 506	Fluran® F-5500-A extension tube	0.51	0.91	10	1
952 20 510	Fluran® F-5500-A extension tube	0.89	0.86	10	1
952 20 512	Fluran® F-5500-A extension tube	1.02	0.86	10	1
952 20 517	Fluran® F-5500-A extension tube	1.42	0.86	10	1
952 20 518	Fluran® F-5500-A extension tube	1.52	0.86	10	1
952 20 522	Fluran® F-5500-A extension tube	2.06	0.86	10	1
952 20 525	Fluran® F-5500-A extension tube	2.79	0.86	10	1

952 22 501	Tygon® LMT-55 extension tube	0.13	0.91	10	1
952 22 502	Tygon® LMT-55 extension tube	0.19	0.91	10	1
952 22 503	Tygon® LMT-55 extension tube	0.25	0.91	10	1
952 22 506	Tygon® LMT-55 extension tube	0.51	0.91	10	1
952 22 509	Tygon® LMT-55 extension tube	0.76	0.91	10	1
952 22 510	Tygon® LMT-55 extension tube	0.89	0.86	10	1
952 22 512	Tygon® LMT-55 extension tube	1.02	0.86	10	1
952 22 515	Tygon® LMT-55 extension tube	1.22	0.86	10	1
952 22 517	Tygon® LMT-55 extension tube	1.42	0.86	10	1
952 22 518	Tygon® LMT-55 extension tube	1.52	0.86	10	1
952 22 522	Tygon® LMT-55 extension tube	2.06	0.86	10	1
952 22 525	Tygon® LMT-55 extension tube	2.79	0.86	10	1
952 22 526	Tygon® LMT-55 extension tube	3.17	0.86	10	1