



Fluovac System and Accessories

- Actively absorbs anesthetic vapors from work area protecting personnel from exposure
- Specially designed nose cones allow optimal access to animals during surgery
- Four different mask sizes available for mice through rabbits (see next page)
- Incorporates special double-mask design to supply and scavenge anesthetic gasses simultaneously
- Requires minimal maintenance
- Supplied with size 1 dual arm (rat)
- Fluosorber canisters ordered separately

Specifications	
Overall Dimensions (H x W x D)	180 x 275 x 300 mm (7.1 x 10.8 x 11.8 in), base unit only
Weight (Including Fluosorber)	8 kg (17.6 lbs)
Weight (Excluding Fluosorber)	7 kg (15.4 lbs)
Power	240 VAC, 50 Hz; 110 VAC, 60 Hz; 37 Watts

Order #	Product
ST1 34-0387	Fluovac, 240 VAC, 50 Hz
ST1 34-0388	Fluovac, 110 VAC, 60 Hz

ST1 34-0415



Fluosorber

- Activated charcoal canister that is packed to optimize halogenated gas absorption
- Stand-alone for passive scavenging or used with the Fluovac for efficient active scavenging
- Canister weights 1,200 g (2.65 lbs) before use and 1,400 g (3.1 lbs) when full
- Canister is designed to scavenge 200 g of anesthetic agent and then can be thrown away in regular trash or incinerated in supplied bags
- Fluosorber is factory tested with an air flow of 60 LPM to ensure resistance is not above 2 cmH₂O

Order #	Product
ST1 34-0415	Fluosorber Canisters, case of 6



ST1 34-0392

Fluovac Exhaust Tube Mount

- 22 mm tapered adapter that fits into the exhaust port on the Fluovac Base Unit in order to adapt to tubing allowing the user to vent gasses like nitrous oxide
- Fits with item ST1 34-0426 corrugated plastic tubing

Order #	Product
ST1 34-0392	Fluovac Exhaust Tube Mount



Fluovac Mask Support

- Used to hold Fluovac Masks in place on a flat surface

Order #	Product
ST1 34-0394	Fluovac, Mask Support for Size 1 Masks
ST1 34-0395	Fluovac, Mask Support for Size 2 Masks
ST1 34-0396	Fluovac, Mask Support for Size 3 Masks