



VWR® Traceable® Precision Lollipop Thermometers

Your lollipop just got sweeter

- Rapid response of 2 seconds within 1 degree, °F/°C switchable
- Water/shockproof IP64 case provides flawless readings in water baths, incubators, tanks
- Stainless-steel probe efficiently monitors temperature in liquids, semi-solids, air/gases
- Supplied: battery, suction cups, magnet, Velcro®, Traceable® Certificate

Applications

Use to take temperature readings, with min and max memory for efficient documentation.

VWR® Traceable® Precision Lollipop Thermometers

Range	Resolution	Accuracy	Sensor Type	Cat. No.	UOM
-58 to 482°F (-50 to 250°C)	0.1°	±0.5°C (-20 to 200°C), ±1°C otherwise	Stainless-steel stem type	620-2722	1 ea.
		±0.5°C (-20 to 150°C), ±1°C otherwise	Stainless-steel probe with cable	620-2723	

Individually calibrated ISO 17025 certified

Individually certified—traceable to NIST