

Product Specification

Meets ACS Specifications.

Test	Specification
Assay	$\geq 97 \%$
Insoluble Matter	$\leq 0.005 \%$
Nitrates (NO ₃)	$\leq 0.001 \%$
pH of 5% Solution at 25°C	6.7–7.3
Residue after Ignition	$\leq 0.01 \%$
Sulphate (SO ₄)	$\leq 0.001 \%$
Trace Impurities (in ppm)	
Chloride (Cl)	≤ 5
Heavy Metals (as Pb)	≤ 5
Iron (Fe)	≤ 5

For Laboratory, Research or Manufacturing use.



Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000