

VWR[®] TraceableGOTM Data Logger Thermometers with external probes





VWR TRACEABLEGO™ BLUETOOTH® DATA
LOGGER THERMOMETERS

- Use free app or with TraceableLIVE®
- Place device (ambient only) or probe inside any container
- 21 CFR 11 compliant with TraceableLIVE® subscription
- External bottle probe: temperature-buffered probe ideal for monitoring conditions of vaccines in transit
- External bullet probe: monitor air and liquids, the sensor and cable may be completely immersed
- External RTD platinum probe: monitor temperatures below –80°C ensuring critical samples are not compromised, 9' FEP cable

TRACEABLE™ CERTIFICATE

An individually-numbered Traceable™ certificate provided with each unit, assures accuracy from an ISO/IEC 17025:2005 (1750.01) calibration laboratory accredited by A2LA. It indicates traceability of measurements to the SI units through NIST or other recognized national measurement institutes (NMI) that are signatories to the CIPM mutual recognition agreement.

Track environmental conditions of reagents, clinical samples, blood, vaccines, tissues, pharmaceuticals, and food items during transport to any location. Bluetooth data loggers monitor shipments while in transit. Compact, durable data loggers fit almost anywhere, even inside a small cooler.

Parameters	Range	Resolution	Accuracy	Probe	Cat. No.
Temperature	–20.00 to 70.00°C (–4.00 to 158.00°F)	0.1°	±0.4°C between –10 and 70°C, otherwise ±0.5°C	"Ambient (no external probe)"	76214-376
	–50.00 to 70.00°C (–58.00 to 158.00°F)	0.1°	±0.3°C	Bullet Probe	76321-888
Humidity & Temperature	0.0 to 95.0% RH, (non-condensing)	0.1°	"±3% RH between 5 to 75%, otherwise ±5% RH"	"Ambient (no external probe)"	76214-378
	–20.0 to 70.0°C (–4.0 to 158.0°F)	0.1°	±0.4°C between –10 and 70°C, otherwise ±0.5°C		
Temperature	–50.00 to 70.00°C (–58.00 to 158.00°F)	0.1°	±0.3°C	Bottle Probe	76321-890
	–90.00 to 100.00°C (–130.00 to 212.00°F)	0.1°	±0.2°C	Platinum Probe	76321-892



Upgrade with TraceableLIVE® subscription

TraceableGO™ App Free Download



FREE

Features

Configure datalogger:

- Set alarm
- Change logging interval
- Toggle between °C/°F
- Start/stop options
- Memory wrap

Export, email and save PDF to mobile device

21 CFR Part 11 compliance

Unlimited cloud data storage

Data download to CSV, or secured PDF

Graphical UI to analyze data

Summary data display: min/max, kinetic mean, time in alarm

Storage of trip parameters

TraceableLIVE® App Free Download



\$30/MONTH

Features

Configure datalogger:

- Set alarm
- Change logging interval
- Toggle between °C/°F
- Start/stop options
- Memory wrap

Export, email and save PDF to mobile device

21 CFR Part 11 compliance

Unlimited cloud data storage

Data download to CSV, or secured PDF

Graphical UI to analyze data

Summary data display: min/max, kinetic mean, time in alarm

Storage of trip parameters

Download FREE TraceableGO™ Mobile APP

Use app to configure data loggers, start/stop options, change memory mode, unit preference, alarm setting, logging interval, and enable/disable alarm. When all the parameters are set, simply press start on the logger to begin monitoring. Once the shipment has been received, double press start/stop on logger to enable Bluetooth and export data to your device. You will get an immediate response to see if your contents have been exposed to conditions outside the set parameters.



Setting science in motion to create a better world

From breakthrough discovery to agile delivery, we offer an extensive portfolio of mission-critical products, services, and solutions. We are a trusted global partner to customers in the life sciences, advanced technologies, and applied materials industries.



VWR, part of Avantor, provides an integrated, seamless purchasing experience optimized for your success.

Learn how the new Avantor is moving forward at avantorsciences.com

©2019 Avantor®, Inc. All rights reserved.



800 932 5000 | VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.