

Our trigger sprayers attach to compatible spray bottles for dispensing bottle contents. They have a squeeze-pump handle and an adjustable nozzle to change from a fine mist to a concentrated stream.

- Use with most trigger heads
- Translucent plastic to monitor chemical levels
- Adjustable nozzle




Qorpak offers a variety of styles and sizes:

Bottle with Trigger Sprayer Attached

Capacity oz (ml)	Neck Finish	Dimensions Od x hgt mm	Bottle Resin	Bottle Style	Dip Tube Length	Case Qty.	Item#	VWR Item#
16 (480)	28-400	70 x 174	Nylon/PE	Modern Round	6 5/8"	6	PLC-06934	11311-616
16 (480)	28-400	61 x 200	HDPE	Cylinder	7 3/8"	6	PLC-03431	16344-101
22 (650)	28-400	84 x 191	HDPE	Graduated Oval	7 3/8"	4	PLC-03499	16344-508
32 (960)	28-400	88 x 198	Nylon/PE	Modern Round	8 1/8"	6	PLC-07581	76487-100
32 (960)	28-400	79 x 238	HDPE	Cylinder	9 1/8"	12	PLC-03438	93001-660

Our Mission is to improve our customers' net income through packaging your business:



Over the last 3 years, our customers' net income increased \$200+ million as a direct result of working with us.

Products Offered

- Jugs
- Narrow and Wide Mouth Bottles
- Dispensing Bottles
- Cans
- Phenolic Caps
- Polypropylene Caps
- Beakers
- Flasks
- Vials
- And more

Markets Served

- Pharmaceutical
- Biotech and Life Sciences
- Agricultural
- Chemical
- Education
- Forensic
- Oil Analysis
- Paints and Coatings
- Water/Wastewater
- Medical
- Petroleum and Fuels