

Material	85908.180
Material description	Triethylamine
Grade	GPR RECTAPUR
CAS Number	121-44-8
Molecular formula	(C ₂ H ₅) ₃ N
Molecular mass	101.19

Characteristics	Specifications
Assay	≥ 99 %
IR Spectrum	Passes test
Refractive index (20°C)	1.399 - 1.401
Water	≤ 0.2 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.