

Certificate of Conformity

VWR Catalogue Number	See List
Description	Electroporation cuvettes
Country of Origin	Manufactured in Taiwan ROC
Date of Issue (yyyy-mm-dd)	2025-07-11

Catalogue Number	Description	Volume [µl]
732-1135	Cuvettes with 1 mm gap size, grey cap	20 - 90
732-1136	Cuvettes with 2 mm gap size, blue cap	40 - 400
732-1137	Cuvettes with 4 mm gap size, yellow cap	80 - 800

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Clear medical grade Polycarbonate.
Sterilization	Product is gamma irradiated with a specified dose range of 25-40kGy.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).Product is manufactured from synthetic

material in a clean room environment.

Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Polycarbonate resin products are manufactured using BPA as a monomeric building block and is consumed during the polymerization process to manufacture polycarbonate. Technically unavoidable trace levels (typically less than 100 milligrams per kilogram polycarbonate or 100 parts per million, ppm) may remain in the polymer matrix of all polycarbonate resins and, therefore, we expect trace levels of BPA to be present. Residual BPA is also fully bound in the polycarbonate resin matrix and would not typically be released into air (it is not volatile at room temperature) unless it is exposed to temperatures where polycarbonate would be in the molten state.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	- 20°C to 60°C
Shelf Life	2 years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	03-2022
2	2.0	Modified Quality system compliance and disclaimer	07-2025