

VWR[®] Specimen Solutions

Securing science, one sample at a time



Sample integrity for critical science

VWR® specimen solutions comprise a comprehensive portfolio of specialized, durable specimen containers designed for clinical, pharmaceutical and biotechnology applications, featuring options from transparent polypropylene to secure ClikSeal™ and high-density polyethylene (PE) containers. This, combined with specialized PE caps and dual-click, 95 kPa-compliant sealing technology, ensures sample integrity and secure transport for diverse laboratory needs. Explore the full product range.

Security

ClikSeal™ technology
Dual-click audible locking
Leak-resistant closures

Versatility

High-clarity polypropylene (PP)
Ultrarugged HDPE
Multiple physical formats

Purity

Sterile options
Nonsterile options
Rigorous life science standards

Several specialized lines.
One trusted source.



**VWR® SPECIMEN CONTAINERS
(PP WITH PE CAPS)**

General-purpose, high-clarity daily
laboratory staples.



VWR® SPECIMEN CONTAINERS (PP)

Classic, all-polypropylene configurations
optimized for chemical compatibility.



**VWR® CLIKSEAL™ SPECIMEN
CONTAINERS**

Advanced thread design engineered for
leakproof sample transport.



**VWR® MULTIPURPOSE
CONTAINERS (HDPE)**

Heavy-duty, large-volume, impact-
resistant bulk storage.



**VWR® SPECIMEN CONTAINERS
WITH DUAL-CLICK LOCK**

Ultimate physical safety with audible/
tactile closure verification.

ORDERING INFORMATION

Container Line	Volume (mL/oz)	Material	Cap Style	Sterility	Qty/Case	Cat. No.
VWR Polypropylene Specimen Containers, Shatter-Resistant						
Specimen Container	148 ml (5 oz.)	PP (Polypropylene)	Cap Not Included (Bulk)	Non-sterile	500	89107-308
Specimen Container	148 ml (5 oz.)	PP (Polypropylene)	Cap Not Included (Bulk)	Sterile	500	89125-840
Specimen Container	237 ml (8 oz.)	PP (Polypropylene)	Cap Not Included (Bulk)	Non-sterile	500	89107-310
Specimen Container	237 ml (8 oz.)	PP (Polypropylene)	Cap Not Included (Bulk)	Sterile	500	89125-842
Specimen Container	237 mL (8 oz.)	PP (Polypropylene)	Container w/Lid	Sterile	100	89125-846
Container Accessory	N/A (Lid Only)	PP (Polypropylene)	Lid Only	Non-sterile	500	89107-312
Container Accessory	N/A (Lid Only)	PP (Polypropylene)	Lid Only	Sterile	500	89125-850
VWR Specimen Containers, Polypropylene with PE Screw Caps: Translucent & Graduated						
Specimen Container	120 ml (4 oz.)	PP (Polypropylene)	Twist Cap (Blue)	Sterile	100	15704-092
Specimen Container	133 ml (4.5 oz.)	PP (Polypropylene)	Twist Cap (White)	Sterile	100	15704-085
Specimen Container	133 ml (4.5 oz.)	PP (Polypropylene)	Twist Cap (White)	Sterile	100	15704-088
Specimen Container	133 ml (4.5 oz.)	PP (Polypropylene)	Screw Cap (White)	Sterile	100	25384-078
Specimen Container	133 ml (4.5 oz.)	PP (Polypropylene)	Screw Cap (Purple)	Sterile	100	25384-144
Specimen Container	133 ml (4.5 oz.)	PP (Polypropylene)	Screw Cap (Purple)	Sterile	100	25384-146
Specimen Container	133 ml (4.5 oz.)	PP (Polypropylene)	Screw Cap (Purple)	Non-sterile	500	25384-148
Specimen Container	150 ml (4.5 oz.)	PP (Polypropylene)	Snap Cap (Clear)	Non-sterile	500	15704-087
VWR KlikSeal Specimen Containers, Sterile Polypropylene						
KlikSeal Container	120 ml (4 oz.)	PP (Polypropylene)	KlikSeal Audible Lock	Sterile	200	15704-008
KlikSeal Container	90 ml (3 oz.)	PP (Polypropylene)	KlikSeal Audible Lock	Sterile	200	15704-002
VWR Multipurpose Containers, HDPE						
Multipurpose Container	120 ml (4 oz.)	HDPE	Snap-on Tight Cap	Non-sterile	250	89009-662
Multipurpose Container	250 ml (8.5 oz.)	HDPE	Snap-on Tight Cap	Non-sterile	250	89009-664
Multipurpose Container	473 ml (16 oz.)	HDPE	Snap-on Tight Cap	Non-sterile	100	89009-666
Multipurpose Container	946 ml (32 oz.)	HDPE	Snap-on Tight Cap	Non-sterile	100	89009-668
Multipurpose Container	1.9 L (64 oz.)	HDPE	Snap-on Tight Cap	Non-sterile	50	89009-670
Multipurpose Container	2.5 L (84.5 oz.)	HDPE	Snap-on Tight Cap	Non-sterile	50	89009-672
Multipurpose Container	4.7 L (1.25 gal.)	HDPE	Snap-on Tight Cap	Non-sterile	25	89009-282
Multipurpose Container	5 L (1.3 gal.)	HDPE	Snap-on Tight Cap	Non-sterile	25	89009-674
VWR Dual Click Specimen Containers, 95 kPa Compliant						
Dual-Click Container	60 ml - Natural Cap (48 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	500	10803-950
Dual-Click Container	60 ml - White Cap (48 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	500	10803-952
Dual-Click Container	60 ml - Blue Cap (48 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	500	10803-956
Dual-Click Container	60 ml - Orange Cap (48 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	500	10803-960
Dual-Click Container	90 ml - Blue Cap (48 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	400	10803-968
Dual-Click Container	90 ml - Natural Cap (48 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	400	10803-962
Dual-Click Container	90 ml - Orange Cap (48 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	400	10803-972
Dual-Click Container	90 ml - White Cap (48 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	400	10803-964
Dual-Click Container	90 ml - Green Cap (53 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	400	10803-978
Dual-Click Container	90 ml - Orange Cap (53 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	400	10803-984
Dual-Click Container	90 ml - Natural Cap (53mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	400	10803-974
Dual-Click Container	90 ml - White Cap (53 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	400	10803-976
Dual-Click Container	120 ml - Orange Cap (53 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	300	10803-996
Dual-Click Container	120 ml - White Cap (53 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	300	10803-988
Dual-Click Container	120 ml - Natural Cap (53 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	300	10803-986
Dual-Click Container	120 ml - Green Cap (53 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	300	10803-990
Dual-Click Container	120 ml - Blue Cap (53 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	300	10803-992
Dual-Click Container	120 ml - Yellow Cap (53 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	300	10803-994
Dual-Click Container	120 ml - White Cap (48mm diam.)	PP (Polypropylene)	Dual-Click Audible	Non-sterile	500	10804-000
Dual-Click Container	120 ml - Blue Cap (48mm diam.)	PP (Polypropylene)	Dual-Click Audible	Non-sterile	500	10804-004
Dual-Click Container	120 ml - Natural Cap (48mm diam.)	PP (Polypropylene)	Dual-Click Audible	Non-sterile	500	10803-998
Dual-Click Container	90 ml - White Cap (48 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Non-sterile	400	10804-012
Dual-Click Container	120 ml - Natural Cap (48mm diam.)	PP (Polypropylene)	Dual-Click Audible	Non-sterile	300	10804-034
Dual-Click Container	120 ml - White Cap (53 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Non-sterile	300	10804-036
Dual-Click Container	120 ml - Yellow Cap (53 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Non-sterile	300	10804-042
Dual-Click Container	90 ml - Yellow Cap (48 mm diam.)	PP (Polypropylene)	Dual-Click (Yellow)	Non-sterile	400	10804-018
Dual-Click Container	90 ml - White Cap (53 mm diam.)	PP (Polypropylene)	Dual-Click (White)	Non-sterile	400	10804-024
Dual-Click Container						
Dual-Click Container	90 ml - Natural Cap (53 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	100 (bulk packed)	10804-050
Dual-Click Container	120 ml - Natural Cap (53 mm diam.)	PP (Polypropylene)	Dual-Click Audible	Sterile	100 (bulk packed)	10804-048

vwr.com
Technical support contact number

Consolidate your inventory. Contact your VWR representative today to request validation samples or a bulk quote.



vwr.com

Any prices, product, and/or services details are current when published and subject to change. Offers subject to product availability. Void where prohibited by law or company policy. Any claims concerning products manufactured by third parties are the sole claims of the manufacturer and not those of Avantor or its affiliates. All products are subject to the terms and conditions of sale of Avantor or its affiliates, including the applicable limitations of liability contained therein. The information herein should not be considered a warranty of any kind. Please carefully review any product information sheets and safety information for the product. Trademarks are owned by Avantor Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](https://www.vwr.com) to view our privacy policy and additional disclaimers. © 2026 Avantor Inc. All rights reserved.