



Certificate of Conformity

VWR CENTRIFUGE TUBES, PS, WITH FLAT CAP

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Catolog No.	Capacity	Version	Cap shape	Tube shape	Packaging configuration
525-0638	15ml	Non Sterile	Flat cap	Conical bottom	BULK 500/CS
525-0639	15ml	Sterile	Flat cap	Conical bottom	25/RACK, 2RACKS/PK, 20 RACKS/CS, 500/CS
525-0640	15ml	Sterile	Flat cap	Conical bottom	50/BAG, 10 BAG/CS 500/CS

Maximum RCF: 3000 g


Materials:

Tube: PS (Polystyrene) Color : Clear

Cap: HDPE (High Density Polyethylene) Color : White

Gasket : PE (Polyethylene) Color : White

Certificate of Conformity



Other Characteristics:

- Made from high quality USP classic VI polystyrene.
- Molded graduations stand out distinctly on these optical clear polystyrene tube.
- Colorless and transparent are perfect for biological reaction experiments.
- Engraved graduation at the body of each tube.
- The flat surface of the polypropylene cap provide an at-a-glance identification area.
- Tubes can be centrifuged to 3000 RCF.
- Not autoclavable.
- Dnase/Rnase free and Non-pyrogenic

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **November 2017**