



Specification

1.00795.2500 n-Hexane for gas chromatography MS SupraSolv®

Specification		
Purity (GC)	≥ 98.0	%
Sum of hexane isomers + methylcyclopentane (GC)	≥ 99.8	%
Identity (IR)	conforms	
Evaporation residue	≤ 3.0	mg/l
Water	≤ 0.01	%
Color	≤ 10	Hazen
GC/MSD (retention range n-undecane to n-tetracontane, scanning area 30 - 600 amu, individual signals (n- tetradecane standard))	≤ 3	ng/ml

Suitable for residue analysis.

Recommended for analysis of polychlorinated Dibenzodioxins and polychlorinated Dibenzofurans (PCDD/F).

Jeannette David

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.