

Material Number: 87041.180
Material description: Ion chromatography cation multi-element standard solution VII
Grade: ARISTAR®
Date specification set: 25th Sep 2019

(All batches of this product after this date are tested to this specification. In the unlikely event that the specification is changed a new certificate will be issued)

Characteristics	Specifications	Traceability
Ba ²⁺	99.8 – 100.2 mg/l	NIST SRM No 3104a
Ca ²⁺	99.8 – 100.2 mg/l	NIST SRM No 3109a
K ⁺	99.8 – 100.2 mg/l	NIST SRM No 3141a
Li ⁺	99.8 – 100.2 mg/l	NIST SRM No 3129a
Mg ²⁺	99.8 – 100.2 mg/l	NIST SRM No 3131a
Mn ²⁺	99.8 – 100.2 mg/l	NIST SRM No 3132
Na ⁺	99.8 – 100.2 mg/l	NIST SRM No 3152a
NH ₄ ⁺	99.8 – 100.2 mg/l	NIST SRM No 194a
Sr ²⁺	99.8 – 100.2 mg/l	NIST SRM No 3153a

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International