



SCIENTIFIC PRODUCTS BEL-ART

Bel-Art Products, Inc.

2024 Broad Street Pocomoke City, MD 21851 USA

Tel: 973-694-0500 * Fax: 973-694-7199

CERTIFICATE OF CONFORMANCE

Product #: 36902-0000 **Mfg. Date:** 2/28/2025 **Lot #:** J392586

Description: "Sterileware" Sterile Sampler Transfer Scoop

Material: High-Impact Polystyrene Natural and White PS Colorant, Sterilized in individual Polyethylene Bags. Materials meet FDA Requirements for Regulation 21 Part 177.1640. White Polystyrene Colorant meets FDA Regulations 21 CFR 177.1640, 178.3297, 175.3.

Other: Irradiation exposure dose (avg.)= 23.5- 50.0 kGy. No testing is done for Pyrogens/Endotoxins at this time. AAMI guidelines are followed. Product should remain in its irradiated state as long as outer packaging remains completely intact. Product should be stored in room temperature dry location.

Based on the information provided to us this product does not contain or come in contact with Lead, Latex, Melamine, Mercury, Formaldehyde, Nitrosamine or other contaminants. To the best of our knowledge based on information provided to us this product does not contain or come in contact with any animal origin and is BSE/TSE protein free. Materials are RoHS, REACH/SVHC, CPSIA Compliant. There is no expiration for this product.

Manufactured by: Bel-Art Products, Inc.

2024 Broad Street Pocomoke City, MD 21851

An ISO-9001:2015 facility (certificate US19/819943409)

Made in the USA

Reviewed and Approved By:

Signature:

Printed Name: Mohamed Mohamed

Title: Quality Assurance Manager

Date: 5/2/2025