

Material	26704.292
Material description	Potassium bromide
Grade	GPR RECTAPUR
CAS Number	7758-02-3
Molecular formula	BrK
Molecular mass	119.00

Characteristics	Specifications
Assay	≥ 98.5 %
Heavy metals (as Pb)	≤ 10 ppm
Cl (Chloride)	≤ 0.5 %
SO ₄ (Sulphate)	≤ 100 ppm
Fe (Iron)	≤ 10 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Jessie Pijck, Head of Laboratory - Haasrode
VWR International bv; Geldenaaksebaan 464; BE-3001
Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.