

Certificate of Conformity

VWR Catalogue Number	See List
Description	Cleanroom Overalls, Sterile
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2021-11-30

Catalogue Number	Colour	Clothing size
113-7881	White	S
113-7882	White	M
113-7883	White	L
113-7884	White	XL
113-7885	White	XXL
113-7886	White	3XL

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

	MD	CD
Break Tensile Strength [N]	> 60	> 35
Elongation rate [%]	> 30	> 30
Trapezoid Testing Strength [N]	> 30	> 30

Product Specifications

The product is in conformance with the requirements of the following standards and regulations:

Material	60g SMS (Spunbonded Polypropylene fabric /Meltblown /Spunbonded Polypropylene fabric).
Sterilization	Products are gamma irradiated according to 11137:2006 and meets a minimum requirement of 10 ⁻⁶ SAL with a dose range of 27.6 – 41.9 kGy.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The above products are not made from natural rubber latex.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Store in room temperature, dry and shady.
Shelf Life	5 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature