



Material	FOR0150AF59001
Material description	Formaldehyde 4 % aqueous solution, buffered
Grade	FCC Formalin 10%
CAS Number	50-00-0
Molecular formula	HCHO
Molecular mass	30.03

Characteristics	Specifications
Appearance	Colourless liquid
Test formaldehyde	Passes test
Density (20/4)	1.012 - 1.014
pH (20°C)	6.8 - 7.2

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Marta Gorewoda - Head of Laboratory - Gliwice
Avantor Performance Materials Poland S.A.; Sowinskiego
11; 44-101 Gliwice

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.