



Certificate of Conformity

VWR CALLIPERS, DIGITAL

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No.	Description
819-0012	VWR Callipers, digital, 0-150 mm
819-0013	VWR Callipers, digital, 0-200 mm

Weight: 255g

Resolution: 0.01 mm

Accuracy: ± 0.03 mm

Other Characteristics:

- Made from stainless steel. Battery life is about 10,000 hours.
- Ergonomic body and friction thumb roller make readings fast.
- Easy to read, jumbo digit display.
- Moisture and dust resistant steel.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **March 2016**