

Certificate Of Quality

Whatman GmbH

Product Information

Product Number: 7195-004
Product Name: WCN WH 5.0 µm
Lot Number: A25229263
Converting Lot: A28732268
Format: Ø 47 mm

Tests

Test	Range
Total thickness [µm]	100 – 140
Air flow rate [s/100ml]	35 – 70
Bubble Point [bar]	0.50 – 0.70

Conformance & Quality systems statement

This is to certify that this product conforms to GE Healthcare Life Sciences specifications.

All products are manufactured via a Quality System certified to **ISO 9001** and **EN ISO 13485** and tested in accordance with documented quality procedures and approved as a result of meeting the required specification.

Quality data are stored.

Electronic signature

This document has been electronically produced and is valid without a signature.

Dassel, 27-Mar-2019 / ABa