

MICROFLEX[®]
ADVANCED DISPOSABLE HAND PROTECTION

93-243

**16 inch / 400 mm Extended Cuff
Disposable Nitrile Gloves**

Extra Long Cuff for Superior Forearm and Product Protection

- Elbow-length cuff (16 inches/400 mm) provides additional protection to the upper arm and extra security against cross contamination
- Made with soft nitrile formulation for comfortable fit and feel while improved cuff design helps prevent gloves from rolling down
- Textured fingers and palm ensure a secure grip in a variety of applications
- Nitrile formulation offers greater puncture resistance and chemical resistance than natural rubber latex
- Gloves are packed in polybags to help protect product integrity



Industries

- Automotive
- Chemical (Basic, Specialty & Petrochemical)
- Electronics (Class 1000 & higher)
- Food Processing
- Life Sciences (Pharmaceutical, Biopharmaceutical & Medical Devices)
- Recycling & Waste Disposal

Applications

- Analytical testing
- Handling compression press and tablet manufacturing
- Handling of paint/glue guns & sprayers
- Industrial cleaning
- Intricate parts handling
- Lab work, blending, compounding, filling
- Maintenance and equipment clean up

TECHNICAL DATA SHEET:

PRODUCT INFORMATION:

	93-243
Material	Nitrile
Color	Blue
Glove Design	Ambidextrous, Powder-Free, Textured Fingers & Palm
Cuff	Beaded Cuff
Manufacturing/QMS Audit Standards	ISO 9001
Regulatory/Standards Compliance	COMPLIANT WITH ASTM D3577, TYPE II, COMPLIANT WITH U.S. FOOD CONTACT REGULATIONS
Packaging	100 gloves per polybag, 10 polybags in a shipper box
Storage	Keep out of direct sunlight; store in a cool and dry place.
Country of Origin	Malaysia

PHYSICAL PROPERTIES:

	TYPICAL VALUES						TESTING METHOD
	5.5 - 6 XS	6.5 - 7 S	7.5 - 8 M	8.5 - 9 L	9.5 - 10 XL	10.5 - 11 XXL	ASTM D3767/EN 420
Length (mm/inches)	400/16	400/16	400/16	400/16	400/16	400/16	
Average Palm Width (mm)	76	86	96	108	120	132	

Freedom from Holes (Inspection level I)	1.5 AQL	ASTM D6319-10/ ASTM D5151-06/EN 374-2
Palm Thickness Single Wall	(mm : 0.12) / (mil : 4.7)	ASTM D3767/EN 420

	BEFORE AGING	AFTER AGING	
Ultimate Tensile Strength (MPa)	> 14	> 14	ASTM D6319-10
Elongation at Break (%)	> 500	> 400	ASTM D6319-10
Force at break (N)	6	6	EN 455-2

ORDERING INFORMATION:

SIZE	5.5 - 6 XS	6.5 - 7 S	7.5 - 8 M	8.5 - 9 L	9.5 - 10 XL	10.5 - 11 XXL
Microflex® 93-243	VWR Catalog Number (US)	75795-894	75795-896	75795-898	75795-900	75795-984
	VWR Catalog Number (CA)	CA75795-894	CA75795-896	CA75795-898	CA75795-900	CA75795-984
	Manufacturer Cat No.	93243060	93243070	93243080	93243090	93243100

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

Performance Standards



Regulatory Compliance



Ansell

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending © 2017 Ansell Limited. All Rights Reserved.

VWR
We Enable Science

1.800.932.5000
vwr.com

Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2017 VWR International, LLC. All rights reserved.