

Combination of epT.I.P.S.® Pipette Tips with Research® plus Pipettes

Eppendorf Research® plus	0.1 – 10 µL dark gray	0.1 – 10 µL M medium gray	0.1 – 20 µL medium gray	0.5 – 20 µL L light gray	0.1 – 20 µL pearl white	0.5 – 100 µL light yellow	2 – 20 µL yellow	2 – 100 µL yellow	2 – 200 µL yellow	20 – 300 µL orange	50 – 1,000 µL blue	50 – 1,250 µL green	500 – 10,000 µL dark green
epT.I.P.S.®	■	■	■	■	□	■	■	■	■	■	■	■	■
epT.I.P.S.® LoRetention	■		■	■	■				■		■		■
ep Dualfilter T.I.P.S.®	■	■		■	■	■	■	■	■	■	■		■
ep Dualfilter	■			■			■	■		■	■		■
T.I.P.S.® LoRetention													
ep Dualfilter T.I.P.S.® SealMax	■			■				■	■	■	■		
Fixed volume, single channel pipettes													
10 µL	■	■	■	■									
20 µL	■	■	■	■									
10 µL	■												
20 µL	■												
25 µL/50 µL/100 µL	■												
200 µL	■												
200/250/500/1,000 µL	■												
Variable volume, single-, 8- and 12-channel pipettes													
0.1 – 2.5 µL	■	■	■										
0.5 – 10 µL	■	■	■	■									
2 – 20 µL	◆	◆	◆	■									
2 – 20 µL								■	■	◆			
10 – 100 µL	■							■	■	■			
20 – 200 µL	■							◆	■	■			
30 – 300 µL	■							◆	◆	■			
100 – 1,000 µL	■										■	■	■
120 – 1,200 µL	■											■	■
0.25 – 2.5 mL	■												
0.5 – 5 mL	■												
1 – 10 mL	■												
Variable volume, 16- and 24-channel pipettes													
1 – 20 µL	□				■								
5 – 100 µL	■												

■ normal volume ◆ limited volume

Combination of epT.I.P.S.® Pipette Tips with Reference® 2 Pipettes

Eppendorf Reference® 2	0.1 – 10 µL S dark gray	0.1 – 10 µL M medium gray	0.1 – 20 µL medium gray	0.5 – 20 µL L light gray	2 – 20 µL yellow	2 – 100 µL yellow	2 – 200 µL yellow	20 – 300 µL orange	50 – 1,000 µL blue	50 – 1,250 µL green	50 – 1,250 µL dark green
epT.I.P.S.®	■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S.® LoRetention	■		■	■				■	■	■	■
ep Dualfilter T.I.P.S.®	■	■		■	■	■	■	■	■		■
ep Dualfilter	■			■		■		■	■		
T.I.P.S.® LoRetention											
ep Dualfilter T.I.P.S.®	■			■		■	■	■	■		
SealMax											
Fixed volume											
1 µL/2 µL	■	■	■								
5 µL/10 µL	■	■	■	■							
20 µL			■	■	■	■	■				
10 µL	■				■	■	■	■			
20 µL	■				■	■	■	■			
25 µL/50 µL/100 µL	■					■	■	■			
200 µL	■						■	■	■	■	■
200/250/500/1,000 µL	■								■		
2 mL/2.5 mL	■										
Adjustable volume											
0.1 – 2.5 µL	■	■	■								
0.5 – 10 µL	■	■	■	■							
2 – 20 µL	◆	◆	◆	■							
2 – 20 µL	■			■	■	■	■	◆			
10 – 100 µL	■				◆	■	■	■			
20 – 200 µL	■				◆	■	■	■			
30 – 300 µL	■				◆	◆	◆	■			
100 – 1,000 µL	■								■	■	■
0.25 – 2.5 mL	■										
0.5 – 5 mL	■										
1 – 10 mL	■										

■ normal volume ◆ limited volume

Combination of epT.I.P.S.® Pipette Tips with Xplorer®/Xplorer® plus Pipettes

Eppendorf Xplorer [®] /plus	0.1 – 10 µL	0.1 – 10 µL M	0.1 – 20 µL	0.5 – 20 µL L	0.1 – 20 µL	0.5 – 100 µL	2 – 20 µL	2 – 100 µL	2 – 200 µL	20 – 300 µL	50 – 1,000 µL	50 – 1,250 µL	50 – 1,250 µL
Eppendorf Xplorer [®] plus	dark gray	medium gray	medium gray	light gray	pearl white	light yellow	yellow	yellow	yellow	orange	blue	green	dark green
epT.I.P.S.®	■	■	■	■	□	■	■	■	■	■	■	■	■
epT.I.P.S.® LoRetention	■		■	■	■	■			■	■			■
ep Dualfilter T.I.P.S.®	■	■		■	■	■	■	■	■	■	■		■
ep Dualfilter	■			■				■		■	■		
T.I.P.S.® LoRetention													
ep Dualfilter T.I.P.S.®	■			■				■	■	■	■		
SealMax													
Variable volume, single-, 8- and 12-channel pipettes													
0.5 – 10 µL	■	■	■	■									
1 – 20 µL	◆	◆	■	■									
5 – 100 µL							◆			◆			
10 – 200 µL							◆	◆	■	■			
15 – 300 µL							◆	◆	■	■			
50 – 1,000 µL							◆	◆		■	■	■	■
50 – 1,200 µL											■	■	■
0.125 – 2.5 mL													
0.2 – 5 mL													
0.5 – 10 mL													
Variable volume, 16- and 24-channel pipettes													
1 – 20 µL	□				■								
5 – 100 µL	■					■							

■ normal volume ◆ limited volume