

BAKERBOND™ XWP 500 Poly CSXB-35

SPE

35µm, 500Å

Hydrophilic Polymer-CSX



Material No.: 7572-01

:

Revision No.: 0

Product Specification

Test	Specification
Appearance	Passes Test
Mean Particle Size µm	30 – 40
Carbon (C), For Information Only	%
Nitrogen (N)	3.5 – 5.5 %
Sulfur (S) (Target 2.5 – 3.0%)	≥ 2.0 %
Pressure (psi)	Passes Test
Liquid Chromatography Use Test	Passes Test
Glucose Oxidase, For Information Only	
IgG, K', For Information Only	
Lysozyme, K', For Information Only	
Lysozyme Retention Time , minutes	16 – 18
Internal Control – g Globulin Capacity (by Frontal Chromatography)(Rabbit or human IgG), For Information Only	
Internal Control – Bacterial Count, For Information Only	CFU/ml

For drug manufacturing or R&D use only. Not on TSCA Inventory.

Storage Condition: Store between 2°C and 15°C

Packaging Site: Phillipsburg Mfg Ctr & DC

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700

Avantor Performance Materials, LLC

100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700

Page 1 of 1