

Isopropyl Alcohol 70% Solution, Sterile
BAKER
Made with USP Water for Injection



Material No.: P002-03
Revision No.: 0

Product Specification

GMP Manufactured Product

Test	Specification
USP – Specific Gravity at 20°/20°C	0.872 – 0.883
Appearance	Passes Test
USP – Acidity (mL)	≤ 1.0
USP – Nonvolatile Residue (mg)	≤ 5.0
USP – Assay	68.0 – 72.0 %
Irradiation: Minimum Delivered Dose (K Gy)	25.0 – 50.0
Irradiation: Maximum Delivered Dose (K Gy)	25.0 – 50.0
Irradiation: Specified Dose Range (Lower)	25.0 – 50.0
Irradiation: Specified Dose Range (Upper)	25.0 – 50.0

GMP Manufactured Product
CAUTION: For Manufacturing, processing or repackaging
Irradiation Date
Irradiation Batch Number
The above referenced lot has passed the current USP Test for Sterility.

Packaging Site: Paris Mfg Ctr & DC

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700

Avantor Performance Materials, LLC

100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700

Page 1 of 1