

Product Specification / Certificate of Compliance

Specification	Product Code	JPB5000P
	Lot Number (if applicable)	M04321
	Description	5000ml BEAKER WITH HANDLE WITH PRINTED GRADUATIONS
	Material	POLYPROPYLENE
	Additional materials (if applicable)	
	Nominal Capacity (ml)	5000
	Pack quantity	

Dimension*	Height (mm)	253
	Width (mm)	218mm INCLUDING SPOUT NOT INCLUDING HANDLE
	Length(mm)	N/A
	Diameter (mm)	200mm NOT INCLUDING SPOUT NOT INCLUDING HANDLE
	Neck Type	N/A

Other	Brimful Capacity (ml)	5800
	Weight (g)	355±15
	Colour	NATURAL

The above product is compliant to the following standards only if marked 'Yes'

Compliance	Food safe (EU/FDA)	YES
	BSE / TSE free	YES
	REACH compliant	YES
	Latex free	YES
	RoHS compliant	UNKNOWN

*All physical dimensions are $\pm 3\%$ and weight is $\pm 1g$, unless otherwise stated.

SciLabware Limited hereby confirms that the enclosed goods have been produced and inspected to all relevant specifications and procedures laid out in our Quality System, in accordance to our registration to BS EN ISO 9001:2008

Signed on behalf of SciLabware Limited by:

Sandra Bell
 Sales & Marketing Assistant



Issue Date	Wednesday, 19 December 2012
Issue Number	1

