







ISOMETRIC VIEW
SCALE 1:1

UNCONTROLLED COPY
10/27/22




- NOTES:
- VIAL IS SUITABLE FOR FREEZE DRYING.
 - VIAL MEETS THE REQUIREMENTS OF THE USP <660> COMPLIANCE.

THIS CONFIDENTIAL AND PROPRIETARY DRAWING AND THE INFORMATION IN IT ARE THE
EXCLUSIVE PROPERTY OF DUREN WHEATON KIMBLE. IT IS PROVIDED WITH THE ASSUMED
AND NOTWITHSTANDING ANY OTHER NOTICES CONTAINED HEREIN, IT SHALL BE REPRODUCED OR DISCLOSED TO
THIRD PARTIES WITHOUT THE PRIOR WRITTEN CONSENT OF DUREN WHEATON KIMBLE. THIS DRAWING IS
UNCONTROLLED AND FOR REFERENCE ONLY UNLESS OTHERWISE STAMPED.

TOLERANCES
UNLESS OTHERWISE
SPECIFIED
X ± 0.00mm
X.X ± 0.01mm
X.XX ± 0.02mm



ALL ITEMS ARE LOW PARTICULATE CLEANED WITH WFI FINAL RINSE, DEPYROGENATION AND AUTOCLAVE



REV	DESCRIPTION OF CHANGE	DRAWN BY	DATE	ECN
-	INITIAL ISSUE	KW	01/21/20	201904041909400

3	230	SEAL MATERIAL: ALUMINUM .008 ±.0008 BUTTON MATERIAL: BLUE POLYPROPYLENE			
2	230	GRAY BUTYL, FM 457/0			
1	210	TYPE I AMBER BOROSILICATE GLASS, 54 EXPANSION AS PER ASTM E438 SEE NOTE 2			
ITEM	QTY.	MATERIAL: DESCRIPTION AND SPECIFICATIONS			
TITLE:					
2ml Amb CPKit,Lyo Stpr,Blu FlpOff Seal					
ITEM NUMBER:		DRAWING NUMBER:		DRAWN BY:	
W990002-CPA		CSD-W990002-CPA		EL	
				DRAWN DATE:	
				04/18/20	
SAP NUMBER:		OLD MATERIAL NUMBER:		SHEET SIZE:	
W990002-CPA		N/A		B	
				SHEET NUMBER:	
				1 OF 1	
				SHEET SCALE:	
				2:1	