

Certificate of Conformity

VWR Catalogue Number	See List
Description	Bottle-Top Vacuum Filtration System, High Performance
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2021-09-14

Complete Filtration Unit, Polyethersulfone Membrane

Catalogue Number	Type	Capacity [L]	Pore Size [µm]	Membrane Diameter Ø [mm]
514-0367	Centrifuge tube filtration unit	15	0.1	40
514-0317	Centrifuge tube filtration unit	15	0.2	40
514-0318	Centrifuge tube filtration unit	15	0.45	40
514-0366	Centrifuge tube filtration unit	50	0.1	40
514-0308	Centrifuge tube filtration unit	50	0.2	40
514-0309	Centrifuge tube filtration unit	50	0.45	40
514-0348	Media bottle filtration unit	250	0.1	66
514-0296	Media bottle filtration unit	250	0.2	66
514-0299	Media bottle filtration unit	250	0.45	66
514-0349	Media bottle filtration unit	500	0.1	90
514-0297	Media bottle filtration unit	500	0.2	90
514-0300	Media bottle filtration unit	500	0.45	90
514-0350	Media bottle filtration unit	1000	0.1	90
514-0298	Media bottle filtration unit	1000	0.2	90
514-0301	Media bottle filtration unit	1000	0.45	90

Filtration Funnels, Polyethersulfone Membrane

Catalogue Number	Type	Capacity [L]	Pore Size [µm]	Membrane Diameter Ø [mm]
514-0363	Filtration funnel for centrifuge tube	15	0.1	40
514-0364	Filtration funnel for centrifuge tube	15	0.2	40
514-0365	Filtration funnel for centrifuge tube	15	0.45	40
514-0359	Filtration funnel for centrifuge tube	50	0.1	40
514-0360	Filtration funnel for centrifuge tube	50	0.2	40
514-0362	Filtration funnel for centrifuge tube	50	0.45	40
514-0354	Filtration funnel for media bottle	250	0.1	66
514-0302	Filtration funnel for media bottle	250	0.2	66
514-0305	Filtration funnel for media bottle	250	0.45	66
514-0355	Filtration funnel for media bottle	500	0.1	90
514-0303	Filtration funnel for media bottle	500	0.2	90
514-0306	Filtration funnel for media bottle	500	0.45	90
514-0358	Filtration funnel for media bottle	1000	0.1	90
514-0304	Filtration funnel for media bottle	1000	0.2	90
514-0307	Filtration funnel for media bottle	1000	0.45	90

Polystyrene Bottles with Polypropylene Caps, Sterile

Catalogue Number	Capacity [L]
391-0286	1000
391-0283	250
391-0284	500

Polycarbonate Bottles with Polypropylene Caps, Sterile

Catalogue Number	Capacity [L]
391-0299	1000
391-0297	250
391-0298	500

Polyethylene terephthalate glycol Bottles with Polypropylene Caps, Sterile

Catalogue Number	Capacity [L]
391-0296	1000
391-0294	250
391-0295	500

Polystyrene Bottles with Polypropylene Caps, Non-Sterile

Catalogue Number	Capacity [L]
391-0317	250

Polycarbonate Bottles with Polypropylene Caps, Non-Sterile

Catalogue Number	Capacity [L]
391-0323	1000
391-0321	250
391-0322	500

Polyethylene terephthalate glycol Bottles with Polypropylene Caps, Non-Sterile

Catalogue Number	Capacity [L]
391-0316	1000
391-0315	500

Quality System Compliance

Products are manufactured under **ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

These Products are manufactured in 100,000 Clean facility that is certified to ISO 13485 Device.

Material	Bottle Polystyrene or Polycarbonate or Polyethylene terephthalate glycol. Cap Polypropylene.
Sterilization	Product are gamma irradiated and meets VD Max SAL 10 ⁻⁶ with a specified dose range of 25-40 kGy (Which is not applicable for Non-Sterile Products).
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Plastic Component materials are non-cytotoxic and meet the Elution Test 10993-5, 1999.
Non-Pyrogenic Statement	The acceptance level for product is less than 20 EU/Device. (LAL Gel Clot Method).
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Phthalates are not deliberately added to the product. This includes di (2-ethylhexyl) phthalate (DEHP), di butyl phthalate (DBP), benzyl butyl phthalate (BBP), diisononyl phthalate (DINP), diisodecyl phthalate (DIDP) and di-n-octyl phthalate (DNOP).
RoHS	Not Applicable.
REACH Statement	Not Applicable.

Storage Conditions Ambient temperature.

Shelf Life 5 years for Sterile Products

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.