

Material	23347.105
Material description	2,6-Dichlorophenolindophenol sodium salt hydrate
Grade	Analytical reagent ACS
CAS Number	1266615-56-8
Molecular formula	$C_{12}H_6Cl_2NNaO_2 \cdot xH_2O$
Molecular mass	290.08

Characteristics	Specifications
IR Spectrum	Passes test
Interfering dyestuffs	Passes test ACS
Loss on drying (120°C)	≤ 12 %
Conforms to ACS	Passes test

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Jessie Pijck, Head of Laboratory - Haasrode
VWR International bv; Geldenaaksebaan 464; BE-3001
Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.