


Terumo Syringe



Product Information Card

Description:

Sterile **3-part syringe** for single use of **Terumo Corporation Japan origin**.

The syringe is intended for the injection of fluids - immediately after filling - into the human body, for aspiration of fluids or for the sampling of blood or other body fluids.

Components:

Does not contain components made of PVC or natural rubber latex.

Description	Raw Materials
Barrel	Polypropylene
Plunger	Polypropylene
Gasket	Thermoplastic elastomer
Lubricant	Silicone oil



Product Range:

Product Code	Size	Tip	replaces
SS*02LE1	2.5 ml	Concentric, Luer Lock	
SS*02SE1	2.5ml	Concentric, Luer Slip	SS+02S1
SS*05LE1	5 ml	Concentric, Luer Lock	
SS*05SE1	5 ml	Concentric, Luer Slip	SS+05S1
SS*10LE1	10 ml	Concentric, Luer Lock	

Technical Data:

Product Code	Nominal Capacity	Graduated Capacity	Increments
SS*02LE1	2.5 ml	2.5 ml	0.1
SS*02SE1	2.5 ml	2.5 ml	0.1
SS*05LE1	5 ml	6 ml	0.2
SS*05SE1	5 ml	6 ml	0.2
SS*10LE1	10 ml	12 ml	0.2

Single unit: blister pack (film/film)

Packaging:

Product Code	Unit Box	Shipping Carton		
		Quantity	Gross Weight	Dimensions (LxWxH)
SS*02LE1	100	2400	6.00 kg	578 x 386 x 436 mm
SS*02SE1	100	2400	6.00 kg	580 x 390 x 440 mm
SS*05LE1	100	1600	6.25 kg	544 x 404 x 444 mm
SS*05SE1	100	1600	6.25 KG	580 X 390 X 440 mm
SS*10LE1	100	800	15.00 kg	580 x 350 x 377 mm

Terumo Syringe

Sterilization:	With e-beam radiation to a Sterility Assurance Level of at least 10^{-6}
Storage:	Shelf life of 5 years <i>* The batch number and expiry date are embossed in the film of the blister pack</i> <i>* Do not store at extreme temperature and humidity</i>
Precautions for Use:	Do not use if package is damaged Do not reuse Dispose of safely after single use to avoid risk of infection
Status of Regulation:	Class IIa, sterile product following MDD 93/42/EEC
Notified body:	TÜV Rheinland 0197
Manufacturer:	TERUMO CORPORATION JAPAN 2-44-1 Hatagaya, Shibuya-ku, Toyko 151-0072 Japan
EC representative:	TERUMO EUROPE Interleuvenlaan 40 3001 Leuven Belgium



Terumo Corporation
+81 3 3374 8111

Terumo Europe NV
+32 16 38 12 11

Terumo Medical Products EMEA
+32 16 38 12 11

EMEA Sales Offices

**Terumo Europe NV
African Business Division**
+32 16 38 13 08

**Terumo Europe NV
Benelux Sales Division
Belgium**
0800 14468
The Netherlands
0800 0231938

**Terumo Europe NV
Emerging Markets Division**
+32 16 38 12 11

Terumo Deutschland GmbH
+49 6196 80 230

Terumo Deutschland GmbH
Switzerland
+41 56 419 10 10

Terumo Europe España SL
+34 902 10 12 98

Terumo France S.A.S.
+33 130 96 13 00

Terumo Italia S.R.L.
+39 6 94 80 28 00

Terumo Russia LLC
+7 495 988 4740

Terumo Sweden AB
+46 3174 85 880

**Terumo Sweden AB
Denmark**
+45 7020 93 80

Terumo Middle East FZE
+971 4 292 0200

Terumo UK Ltd
+44 1276 480 440

**Terumo BCT Tibbi Cihazlar Dagıbm
ve Hizmetleri A.S.**
+90 216 645 92 00

