


Physics

From classic and cutting-edge apparatus and demos, to time-saving activities, durable equipment, and essential supplies for your lab - we have what you need to turn physics lessons into physics connections.



Specifications:

P1035 Multimeter

DC voltage:	200 mV – 600 V, 5 ranges, $\pm 0.5\% \pm 2$ digits
AC voltage:	200/ 600 V, 2 ranges, $\pm 1.2\% \pm 10$ digits
DC current:	2000 μ A – 10 A, 4 ranges, $\pm 1\% \pm 2$ digits
Resistance:	200 Ω – 2000 k Ω , 5 ranges, $\pm 0.8\% \pm 2$ digits
Display:	3½ digit LCD, 27 mm, max: 1999
Operating voltage:	9 V battery
Safety classification:	CAT III 600 V (IEC-1010-1)
Fuse:	F1: F 200 mA / 600 V F2: F 10 A / 600 V, $I_{max.} = 10$ A for 30 s with minimum interval 15 mins
Dimensions:	approx. 70x150x48 mm ³
Weight:	approx. 260 g