



Certificate of Conformity

VWR SINGLE CHANNEL PIPETTES, MECHANICAL, VARIABLE VOLUME, FULLY AUTOCLAVABLE

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat. No.	Capacity(μl)	Accuracy (%)	Imprecision (%)	Pk
613-2891	0.2 - 2	±12.0 - ±1.5	≤6.0 - ≤0.7	1
613-2892	0.5 - 10	±4.0 - ±0.5	≤2.8 - ≤0.4	1
613-2893	2 - 20	±3.0 - ±0.8	≤1.5 - ≤0.3	1
613-2894	10 - 100	±1.6 - ±0.8	≤0.8 - ≤0.2	1
613-2895	20 - 200	±1.2 - ±0.6	≤0.6 - ≤0.2	1
613-2896	100 - 1000	±1.6 - ±0.6	≤0.4 - ≤0.15	1
613-2897	1000 - 5000	±0.6 - ±0.5	≤1.2 - ≤0.15	1
613-2898	1000 - 10000	±2.5 - ±0.5	≤0.6 - ≤0.2	1

Other Characteristics:

- Easily adjustable ejector height system to accommodate virtually all brands of tips.
- UV resistant and fully autoclavable.
- Easy and quick recalibration system.
- Colour-coded shaft for easy volume identification.
- 5000 and 10000 μl version with shaft protected by filter.

The products pass stringent validation & control procedures based on EN ISO 8655. They are calibrated & tested using gravimetry method with distilled water.

The conditions of measurement : Basis of adjustment: Ex; reference temperature: 20°C; relative air humidity: 50%; Barometric pressure:101kpa.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **June 2016**