



ZPACK13950512 PLASTIC DRUM LIQ 200 L UN



SPECIFICATIONS

UN-certification :	1H1/X1.9/200/
Material :	HDPE blue
Filling Opening :	2,000 "
Top Diameter :	581,00 mm
Bottom Diameter :	581,00 mm
Height :	935,000 mm

AVAILABLE ACCESSORIES

Key for opening :	29548.061
Plastic taps :	29548.097
Pump PP Bürkle :	223-1262

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343