


Helping to make recycling simple

Recycle your ART, ART Barrier and SoftFit-L pipette tips, tip boxes and associated plastic packaging

Our pipette tip and box mail-back recycling program provides you with a simple and effective solution for recycling Thermo Scientific™ ART™, Thermo Scientific™ ART™ Barrier and Thermo Scientific™ SoftFit-L™ pipette tips. Simply collect your used pipette tips and boxes and our convenient mail-back program helps take care of the rest.



Recyclable tips

ART barrier hinge rack, sterile packaging option


Tip image	Volume	Description	Tip length (mm)	Tip O.D. (mm)	Tip I.D. (mm)	PK	CS	Sterility	Standard Tip VWR Cat. No.	Low-retention Tip VWR Cat. No.
xS10	10 µL	Standard length, graduated at 2.5 µL, micropoint, universal fit	31.27	5.87	0.406	960	4,800	Sterile	76809-072	76809-070
xL10	10 µL	Extended length, graduated at 2.5 µL, micropoint, universal fit	43.28	5.79	0.432	960	4,800		76809-076	76809-074
sStd20	20 µL	Standard length, non-graduated, micropoint, universal fit	50.42	7.49	0.559	960	4,800		76809-080	76809-078
sStd100	100 µL	Standard length, non-graduated, micropoint, universal fit	50.42	7.49	0.559	960	4,800		76809-060	76809-058
sStd200	200 µL	Standard length; graduated at 50, 100 and 200 µL; micropoint; universal fit	58.83	7.04	0.533	960	4,800		76809-064	76809-062
xL1000	1,000 µL	Extended length; graduated at 100, 200, 500 and 1,000 µL; micropoint; universal fit	88.90	8.89	0.813	768	3,072		76809-084	76809-082
uXL1000	1,000 µL	Extended length, non-graduated, micropoint, universal fit	118.49	8.66	0.711	800	3,200		76809-068	76809-066

For high-value samples where contamination is a concern, ART Barrier pipette tips enable contamination protection and help eliminate costly downtime for decontamination efforts. Ideal for PCR, qPCR, DNA, RNA and any other assays where reliable and repeatable results are required. Cross-contamination can lead to reporting errors and cause irreparable damage to a lab's credibility. The ART Barrier filter is instrumental in helping to eliminate contamination even in the rare event of erroneous excess-volume pipetting.

ART non-barrier rack, sterile packaging option

Tip image	Volume	Description	Tip length (mm)	Tip O.D. (mm)	Tip I.D. (mm)	PK	CS	Sterility	Standard Tip VWR Cat. No.	Low-retention Tip VWR Cat. No.
xS10	10 µL	Standard length, graduated at 2.5 µL, micropoint, universal fit	31.27	5.87	0.406	960	4,800	Sterile	76809-102	76809-100
xL10	10 µL	Extended length, graduated at 2.5 µL, micropoint, universal fit	43.28	5.79	0.432	960	4,800		76809-106	76809-104
ex20	20 µL	Extended length, non-graduated, micropoint, 10 µL tip collar design, Eppendorf™ fit	45.75	6.83	0.457	960	4,800		76809-110	76809-108
sStd200	200 µL	Standard length, non-graduated, micropoint, universal fit	50.42	7.49	0.559	960	4,800		76809-114	76809-112
xL1000	1,000 µL	Extended length; graduated at 100, 200, 500 and 1,000 µL; micropoint; universal fit	88.90	8.89	0.813	768	3,072		76809-098	76809-096

ART non-barrier pipette tips are ideal for everyday lab applications where contamination is not a concern, but accuracy and precision are still required for reliable results.



Our pipette tip and box mail-back recycling program provides you with a simple and effective solution for recycling. Simply collect your used pipette tips and boxes, and our convenient mail-back program helps take care of the rest.

SoftFit-L recyclable tips			
Volume	Description	Packaging	VWR Cat. No.
20 µL	SoftFit-L 20 µL sterile, filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	76809-088
200 µL	SoftFit-L 200 µL sterile, filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	76809-090
300 µL	SoftFit-L 300 µL sterile, filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	76809-086
1,000 µL	SoftFit-L 1,000 µL sterile, filtered tips	96 tips/rack, 8 racks/pk, 4 pk/cs	76809-092
1,200 µL	SoftFit-L 1,200 µL sterile, filtered tips	96 tips/rack, 8 racks/pk, 4 pk/cs	76809-094
20 µL	SoftFit-L 20 µL sterile, non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	76809-116
200 µL	SoftFit-L 200 µL sterile, non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	76809-126
300 µL	SoftFit-L 300 µL sterile, non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	76809-120
1,000 µL	SoftFit-L 1,000 µL sterile, non-filtered tips	96 tips/rack, 8 racks/pk, 4 pk/cs	76809-130
1,200 µL	SoftFit-L 1,200 µL sterile, non-filtered tips	96 tips/rack, 8 racks/pk, 4 pk/cs	76809-094
20 µL	SoftFit-L 20 µL non-sterile, non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	76809-118
200 µL	SoftFit-L 200 µL non-sterile, non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	76809-128
300 µL	SoftFit-L 300 µL non-sterile, non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	76809-122
1,000 µL	SoftFit-L 1,000 µL non-sterile, non-filtered tips	96 tips/rack, 8 racks/pk, 4 pk/cs	76809-132
1,200 µL	SoftFit-L 1,200 µL non-sterile, non-filtered tips	96 tips/rack, 8 racks/pk, 4 pk/cs	76809-124

FedEx is a registered trademark of Federal Express Corporation. Eppendorf is a registered trademark of Eppendorf AG.