



Saniclave-Autoclave 102 - FDA Cleared



Great For: General Medical Dental Chiropody Beauty Body Art Veterinary General lab use Academia

Specifications

Pressure: 15 psi

Internal Temperature: 121°C

Sterilization Cycle: 30 minutes

Capacity: 8L

Net Weight: 5.4kg (12lbs)

Power Requirements: 110-130 VAC, 50-60 Hz, 1250 W

Display: Digital time and/or temperature

Durability: Rugged stainless steel chamber and cover

Safety: Over temperature alarm

Load type: Solid, nonlumened, unwrapped medical instruments

Dimensions External: 13" tall x 12.5" wide x 15" deep Internal: 9" diameter x 8" tall

Weight	<i>18 lbs</i>
Dimensions	<i>17 × 17 × 17 in</i>