



Certificate of Conformity

VWR ALARM THERMOMETER WITH TIMER TRACEABLE

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No.	Range (°C)	Probe	Pk
620-1662	0....200 °C	Curved Stainless Steel	1

Curved Stainless Steel Probe : (L×Ø x W: 70×76x12.7 mm) with three metre cable and internal sensor

Accuracy (°C): ± 2

Resolution (°C): 1

Dimensions: 2 ¾ x 3x ½ inches

Weight: 4.2 ounces

Other Characteristics:

This thermometer, with timer, is designed to monitor liquids, air, gas or semi-solids.

The metal cable may be placed in ovens and incubators making this unit ideal for use in these environments.

Display shows time, current temperature and alarm setting.

Adjustable display and large digits

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **August 2018**