

Material No.: 76163-604

Material Description: VWR syringe tips, classic style, 1.25 mL, sterile, 100/pack

Lot No.: 1005351

Expiration Date: 11/30/2026

Compliance and Procedure Result

Quality System Compliance

Products are manufactured in ISO 13485:2016 and ISO 9001:2015 certified plant. Products are inspected and controlled in accordance with the current applicable product specification.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material Tips: Polypropylene

Manufacturing Environment: Robotic tips are manufactured under clean room conditions according to DIN EN ISO 14644-1:2016-06.

Sterilization: Products are EB (Electron beam) irradiated and dose released upon ISO 11137-2:2013 recommended practices in effect at the time of validation meet the requirement of 10⁶ SAL.

ATP Assay: Tested with detection level < 1*10⁻¹²

PCR inhibitors Free: Product yield equivalent to control < 10 amplifiable targets

Human DNA Free: Tested with detection level < 1 pg

DNase & RNase Free: Not applicable

BSE/TSE: No use of any raw materials produced from or substances derived from animal origin. The manufacturing process of the product does not intentionally use any ingredient of animal origin and no material derived from or exposed to animals affected by.

Country of Origin: Germany

Signed on behalf of VWR By: James Osborn
Quality Employee

5/9/2023