

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Tissue Grinders, Potter-Elvehjem Type
Country of Origin	Manufactured in United Kingdom
Date of Issue (yyyy-mm-dd)	2023-08-14

Glass Vessel

Catalogue Number	Capacity [mL]	Ø ext x Length [mm]	Total Length [mm]
432-0200	2	8 × 120	120
432-0201	5	12 × 135	135
432-0202	10	15 × 150	150
432-0203	15	19 × 155	155
432-0204	30	25 × 175	175
432-0205	50	32 × 195	195

Plain Plunger

Catalogue Number	Capacity [mL]	Total Length [mm]
432-0206	2	230
432-0207	5	235
432-0208	10	270
432-0209	15	270
432-0210	30	270
432-0211	50	270

Plunger with serrated tip

Catalogue Number	Capacity [mL]	Total Length [mm]
432-0212	2	230
432-0213	5	235
432-0214	10	270
432-0215	15	270
432-0216	30	270
432-0217	50	270

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Virgin Polytetrafluoroethylene (PTFE) Grade TF1620. Raw Material complies with Directive FDA 21 CFR 177.1550 and complies with USP Class VI - USP 27, NF 22,2004 for plastic at 70 °C. The test is carried out by 3M-Dyneon who manufactures PTFE used and refers to PTFE raw material. 3M™ DYNEON™ TF1620 – COMPLIANCE STATUS. Glass Homogeniser: Borosilicate Glass Plain Plunger: TF1620 Plunger with serrated tip: TF1620 and Type 304 Stainless Steel
Sterilization	Non-Sterile.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2023-06-14 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Clean and Dry.
Shelf Life	Indefinite when stored under ideal conditions.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	08-2023