



VWR® Single Well Tissue Culture Plates

Description: Single Well, Surface Treated/Non-treated, Sterilized

Purpose: For use in applications where cell attachment is desired or not desired

Materials


Dish: GPPS (General polystyrene) Colour: Clear

Lid: GPPS (General polystyrene) Colour: Clear

Features


- Single well plates have increased internal plating dimension, providing larger area for cell growth
- Plates have a standard external footprint, making them compatible with other instruments that accept SBS-standard plates
- Suitable for robotic handling
- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- DNase/RNase free, Non-pyrogenic
- Autoclavability: No
- Assembled plate and lid height: 17.5 mm
- Max. volume (single well): 90.0 ml
- Working volume (single well): 35.0 ml
- Growth area (single well): 97.0 cm²
- Every inner bag is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system

Dimensions (Unit:mm)



VWR® Single Well Tissue Culture Plates

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Pk	Qty/Case
734-2977	75780-348	Single well Tissue Culture plates, General, Non-treated, Sterile	Yes	10	100
734-2966	75999-254	Single well Tissue Culture plates, Standard, Surface Treated, Sterile	Yes	10	100



VWR® 4 Well Tissue Culture Plates

Description: 4 Well, Surface Treated/Non-treated, Sterilized


Purpose: For use in applications where cell attachment is desired or not desired

Materials

Dish: GPPS (General polystyrene) Colour: Clear

Lid: GPPS (General polystyrene) Colour: Clear

Dimensions (Unit:mm)



Features

- 4-well plates have a pinchbar design, and flat-bottom round wells, with lid
- Lid with notched corners for good fit and easy orientation
- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- DNase/RNase free, Non-pyrogenic
- Available in individual recyclable PET
- Autoclavability: No
- Assembled plate and lid height: 13.5 mm
- Max. volume (single well): 1.86 ml
- Working volume (single well): 1.0 ml
- Growth area (single well): 1.96 cm²
- Every inner bag is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® 4 Well Tissue Culture Plates

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
732-2904	10813-630	4 well Tissue Culture plates, General, Non-treated, Sterile	Yes	1	50	100
732-2905	10814-222	4 well Tissue Culture plates, Standard, Surface Treated, Sterile	Yes	1	50	100



VWR® 6 Well Tissue Culture Plates

Description: 6 Well, Surface Treated/Non-treated, Sterilized


Purpose: For use in applications where cell attachment is desired or not desired

Materials

Dish: GPPS (General polystyrene) Colour: Clear

Lid: GPPS (General polystyrene) Colour: Clear

Dimensions (Unit:mm)



Features


- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- DNase/RNase free, Non-pyrogenic
- Available in individual recyclable PET
- Autoclavability: No
- Assembled plate and lid height: 22.2 mm
- Max. volume (single well): 17.0 ml
- Working volume (single well): 1.9~2.9 ml
- Growth area (single well): 9.6 cm²
- Every inner bag is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



STERILE R

VWR® 6 Well Tissue Culture Plates

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
734-2777	10861-554	6 well Tissue Culture plates, General, Non-treated, Sterile	Yes	1	50	100
734-2323	10062-892	6 well Tissue Culture plates, Standard, Surface Treated, Sterile	Yes	1	50	100



VWR® 12 Well Tissue Culture Plates

Description: 12 Well, Surface Treated/Non-treated, Sterilized


Purpose: For use in applications where cell attachment is desired or not desired

Materials

Dish: GPPS (General polystyrene) Colour: Clear

Lid: GPPS (General polystyrene) Colour: Clear

Dimensions (Unit:mm)



Features


- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- DNase/RNase free, Non-pyrogenic
- Available in individual recyclable PET
- Autoclavability: No
- Assembled plate and lid height: 22.7 mm
- Max. volume (single well): 6.8 ml
- Working volume (single well): 0.76~1.14 ml
- Growth area (single well): 3.85 cm²
- Every inner bag is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



STERILE R

VWR® 12 Well Tissue Culture Plates

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
734-2778	10861-556	12 well Tissue Culture plates, General, Non-treated, Sterile	Yes	1	50	100
734-2324	10062-894	12 well Tissue Culture plates, Standard, Surface Treated, Sterile	Yes	1	50	100



VWR® 24 Well Tissue Culture Plates

Description: 24 Well, Surface Treated/Non-treated, Sterilized


Purpose: For use in applications where cell attachment is desired or not desired

Materials

Dish: GPPS (General polystyrene) Colour: Clear

Lid: GPPS (General polystyrene) Colour: Clear

Dimensions (Unit:mm)



Features


- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- DNase/RNase free, Non-pyrogenic
- Available in individual recyclable PET
- Autoclavability: No
- Assembled plate and lid height: 22.7 mm
- Max. volume (single well): 3.5 ml
- Working volume (single well): 0.38~0.57 ml
- Growth area (single well): 1.93 cm²
- Every inner bag is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



STERILE R

VWR® 24 Well Tissue Culture Plates

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
734-2779	10861-558	24 well Tissue Culture plates, General, Non-treated, Sterile	Yes	1	50	100
734-2325	10062-896	24 well Tissue Culture plates, Standard, Surface Treated, Sterile	Yes	1	50	100



VWR® 48 Well Tissue Culture Plates

Description: 48 Well, Surface Treated/Non-treated, Sterilized


Purpose: For use in applications where cell attachment is desired or not desired

Materials

Dish: GPPS (General polystyrene) Colour: Clear

Lid: GPPS (General polystyrene) Colour: Clear

Dimensions (Unit:mm)




Features

- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- DNase/RNase free, Non-pyrogenic
- Available in individual recyclable PET
- Autoclavability: No
- Assembled plate and lid height: 22.7 mm
- Max. volume (single well): 1.55 ml
- Working volume (single well): 0.19~0.29 ml
- Growth area (single well): 0.83 cm²
- Every inner bag is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® 48 Well Tissue Culture Plates

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
734-2780	10861-560	48 well Tissue Culture plates, General, Non-treated, Sterile	Yes	1	50	100
734-2326	10062-898	48 well Tissue Culture plates, Standard, Surface Treated, Sterile	Yes	1	50	100



VWR® 96 Well Flat Bottom Tissue Culture Plates

Description: 96 Well, Flat Bottom, Surface Treated/Non-treated, Sterilized


Purpose: For use in applications where cell attachment is desired or not desired

Materials

Dish: GPPS (General polystyrene) Colour: Clear

Lid: GPPS (General polystyrene) Colour: Clear

Dimensions (Unit:mm)




Features

- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- DNase/RNase free, Non-pyrogenic
- Available in individual recyclable PET
- Autoclavability: No
- Assembled plate and lid height: 16.0 mm
- Max. volume (single well): 0.39 ml
- Working volume (single well): 0.075~0.2 ml
- Growth area (single well): 0.33 cm²
- Every inner bag is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® 96 Well Flat Bottom Tissue Culture Plates

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
734-2781	10861-562	96 well Tissue Culture plates, Flat Bottom, General, Non-treated, Sterile	Yes	1	50	100
734-2327	10062-900	96 well Tissue Culture plates, Flat Bottom, Standard, Surface Treated, Sterile	Yes	1	50	100



VWR® 96 Well U-Bottom Tissue Culture Plates

Description: 96 Well, U-Bottom, Surface Treated/Non-treated, Sterilized


Purpose: For use in applications where cell attachment is desired or not desired

Materials

Dish: GPPS (General polystyrene) Colour: Clear

Lid: GPPS (General polystyrene) Colour: Clear

Dimensions (Unit:mm)




Features

- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- DNase/RNase free, Non-pyrogenic
- Available in individual recyclable PET
- Autoclavability: No
- Assembled plate and lid height: 16.2 mm
- Max. volume (single well): 0.33 ml
- Working volume (single well): 0.075~0.2 ml
- Growth area (single well): 0.58 cm²
- Every inner bag is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® 96 Well U-Bottom Tissue Culture Plates

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
734-2782	10861-564	96 well Tissue Culture plates, U-Bottom, General, Non-treated, Sterile	Yes	1	50	100
734-2328	10062-902	96 well Tissue Culture plates, U-Bottom, Standard, Surface Treated, Sterile	Yes	1	50	100



VWR® 384 Well Tissue Culture Plates

Description: 384 Well, Surface Treated/Non-treated, Sterilized


Purpose: For use in applications where cell attachment is desired or not desired

Materials

Dish: GPPS (General polystyrene) Colour: Clear

Lid: GPPS (General polystyrene) Colour: Clear

Dimensions (Unit:mm)



Features

- 384-well plates have square wells and with a flat bottom
- Plates have a pinchbar design, and lid with notched corners for good fit and easy orientation
- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- DNase/RNase free, Non-pyrogenic
- Available in individual recyclable PET
- Autoclavability: No
- Assembled plate and lid height: 17.3 mm
- Max. volume (single well): 0.145 ml
- Working volume (single well): 0.01~0.1 ml
- Growth area (single well): 0.1135 cm²
- Every inner bag is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® 384 Well Tissue Culture Plates

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
732-2906	10814-224	384 well Tissue Culture plates, General, Non-treated, Sterile	Yes	1	50	100
732-2907	10814-226	384 well Tissue Culture plates, Standard, Surface Treated, Sterile	Yes	1	50	100