


Fine mesh basket

P/N: 7501135 | SiKo/KT 274/174/36 BR oG DM mD

HUPFER
we make work flow

Technical data



Frame wire thickness:	3.5 mm
Lock wire thickness:	2.0 mm
Mesh size:	1.6 × 1.6 mm
Weight:	0.462 kg
Width:	274 mm
Depth:	174 mm
Height:	36 mm

Similar to illustration, technical modifications reserved. Without decoration.

Rinsing basket with fixed lid, made of high-quality stainless steel wire for storing instruments, lenses and accessories during the cleaning process.

Rinsing basket and lid made of finely-woven stainless steel gauze with frame made of stainless steel wire and secure quick fastener. Base, lid and walls made of wire gauze.

Time and date of the request:
21.03.2023, 10:49:06

All information / dimensions are approximate, technical changes reserved. © Hupfer